



SuperLock Microcentrifuge Tubes



Product Highlights

- Locking tubes for the highest security
- Pierceable cap and frosted writing area
- Rated at 17,000 g
- Certified RNase and DNase free

Product Applications

Sample preparation
Centrifugation
Sample boiling or freezing
Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	See variations
Sterile	No
Autoclavable	Yes
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
RNase free	Yes
Color	Natural
Max. temp.	121 °C
Temperature range	-80 to 121 °C
Temperature min. (SL)	-80 °C
g-Force max.	17,000 x g

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de







More informations about SuperLock Microcentrifuge Tubes

This range of [microcentrifuge tubes](#) offers the ultimate in sample security.

- > Locking cap
- > Piercable cap
- > Graduated
- > Max. centrifugation rate of 0.5 ml and 1.5 ml tubes: 17,000 x g
- > Max. centrifugation rate of 2.0 ml tubes: 13,000 x g
- > Certified RNase and DNase free



All Variations


PRODUCT NAME	PACKAGING SIZE	ART. NO.
 0.5 ml SuperLock Microcentrifuge Tube, Natural Volume: 0.5 ml  Create PDF Data Sheet	1,000 Tubes (1 Bag × 1000 Tubes)	E1405-5100
 1.5 ml SuperLock Microcentrifuge Tube, Natural Volume: 1.5 ml  Create PDF Data Sheet	500 Tubes (1 Bag × 500 Tubes)	I1415-5100

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Cap-Locks for 1.5 / 2.0 ml Microcentrifuge Tubes	100 Pcs. (1 Bag × 100 Pcs.)	I1415-1508

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de