



## Culture Tubes



### Product Highlights

- 14 ml capacity
- White writing area
- Graduations marked at 0.25 ml increments
- Centrifugation rates: 3,000 x g (PP), 1,400 x g (PS)
- Temperature use: -190 °C to 120 °C (PP) / -10 °C to 70 °C (PS)
- Sterile
- Dimensions (HxØ) 95 x 17 mm

## Product Applications

Cell culture

Routine laboratory procedures like centrifugation, heating and freezing

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	See variations
Sterile	Yes
Volume	14 ml
Pyrogen free	Yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about Culture Tubes

These disposable sterile tubes can be used for most routine laboratory procedures. They are biologically inert and exempt from mold release agents. Precision molding with virgin thermoplastics ensures that our tubes will be uniform in size and shape. High resistance to breakage reduces danger in handling infectious or other potentially harmful cultures. Dual position natural caps, which can be left loose for aerobic work or sealed for anaerobic cultures.

Transparent polystyrene tubes will withstand moderate centrifugation speeds (1,400 x g) and temperatures to 70 °C. Translucent polypropylene tubes can be centrifuged at higher speeds (3,000 x g) and resist temperatures from -190 °C to 120 °C.



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">14 ml Polypropylene Culture Tube with Cap</a> <b>Material:</b> Polypropylene  <a href="#">Create PDF Data Sheet</a>	500 Tubes (20 Bags × 25 Tubes)	I1485-0810
 <a href="#">14 ml Polystyrene Culture Tube with Cap</a> <b>Material:</b> Polystyrene  <a href="#">Create PDF Data Sheet</a>	500 Tubes (20 Bags × 25 Tubes)	I1485-2810

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)