



100 µl Bevelled TipOne® Filter Tip, Natural, Refill (sterile), Case



Product Highlights

- Highly effective HDPE filter that contains no additive to interfere with your samples
- Universal tip fit.
- High-quality virgin polypropylene tips, with extremely smooth surfaces for reduced sample retention
- Sterile and Certified RNase, DNase, DNA and Pyrogen free
- Refills dramatically reduce plastic waste compared to pre-racked tips
- Refill top and bottom shells made from 100 % recycled plastic
- All components, excluding film wrap and tips, are recyclable
- Buy by the case to save on packaging and transportation. Reuse the case box to send back your recycling.

Product Applications

DNA applications (e.g. PCR)

RNA applications (e.g. Gene expression analysis)

Protein applications (e.g. Antibody Research)

Cell Culture applications (e.g. Media)

Applications with radio actives

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	S1120-1740-C
Sterile	Yes
Autoclavable	No
Filter tip	Filter
Volume range	1.0 - 100 µl
Volume	100 µl
Pack Size	7,680 Tips (80 Refills × 96 Tips)
Color Code Wafer	Yellow
Tip Type	bevelled
Material	Polypropylene
Graduations	No
Color	Natural
Sustainability attributes	Refills dramatically reduce plastic waste compared to pre-racked tips. Refill top and bottom shells made from 100 % recycled plastic. All components, excluding tips and shrink wrap wrap, are recyclable. Buy by the case to save on packaging and transportation. Reuse the case box to send back your recycling.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



For use with brand + product (SL)

Abimed: Discovery Comfort, LABMATE Optima

Biohit: eLine, ePet, mLine, Picus, Proline

BRAND: Transferpette, Transferpette Electronic, Transferpette Pro, Transferpette S

Capp: CappAero Comfort, CappAero Multi, CappAero Universal, Cappronic Multi, Cappronic Single, ecopipette, Singlechannel

Eppendorf: Reference, Reference 2, Research, Research 3 neo, Research plus, Research pro, Xplorer

Gilson: Pipetman L, Pipetman Neo, Pipetman P, Pipetman Ultra

Hamilton: SoftGrip Pipette, SofTouch

HTL: Discovery Comfort, Discovery Pro, Labmate Pro, Optipette

Nichiryo: 00-NP7, 50F, 50V, NLE, NPLO, NPX

Rainin: Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus

Sartorius: eLine, ePet, mLine, Picus, Picus NxT, Proline, Proline Plus, Tacta

Socorex: Acura 825, Acura 855, Acura electro 925, Acura electro 955, Calibra digital

STARLAB: ErgoOne, ErgoOne E, MicroOne, StarPet

Thermo Fisher Scientific Inc.: Finn Focus Micro, Finn Focus Univ., Finnpipette Digital, Finnpipette Novus Turquoise

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 100 µl Bevelled TipOne® Filter Tip, Natural, Refill (sterile), Case

Nothing goes through!

Protect your samples from unwanted contamination with TipOne® [filter tips](#). Our high-quality filters prevent aerosol cross-contamination without trapping your valuable samples. Plus, they contain no additives that could interfere with your samples.

You can't always see fine aerosols entering your [pipette](#), but they can contaminate other samples later. TipOne® [filter tips](#) block over 99% of aerosols*, making them one of the best filters available.

Sterile TipOne® tips are validated according to EN ISO 11137-2:2007.

The filter tip refill, is an economical, plastic-saving alternative to racked [filter tips](#).

Reduce single-use plastic!

Compared to pre-racked tips, using TipOne® sterile filter refills dramatically reduces the amount of single-use plastic going through your laboratory. The sterile refill is quick and easy to use, just unwrap and place the entire refill into an empty TipOne® rack. The fully enclosed refill means sterility is never compromised. Tip info is printed on every wafer, for full traceability when racked.








Join our journey to sustainable pipetting

TipOne® tip racks are built to last and can be reused again and again with these cost-effective refills. The refill's top and bottom shells are made from 100% recycled plastic. All components of TipOne® sterile filter tips refills are recyclable (excluding tips and shrink wrap). Sign up for our TipOne® recycling service!

*See TipOne® brochure in downloads section for details.



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 100 µl Bevelled TipOne® Filter Tip, Natural, Refill (sterile) Pack Size: 960 Tips (10 Refills × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Refills × 96 Tips)	S1120-1740
 100 µl Bevelled TipOne® Filter Tip, Natural, Refill (sterile), Case Pack Size: 7,680 Tips (80 Refills × 96 Tips)	7,680 Tips (80 Refills × 96 Tips)	S1120-1740-C
 100 µl Bevelled TipOne® Filter Tip, Natural, Racked (sterile) Pack Size: 960 Tips (10 Racks × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	S1120-1840
 100 µl Bevelled TipOne® Filter Tip, Natural, Racked (sterile), Case Pack Size: 7,680 Tips (80 Racks × 96 Tips)  Create PDF Data Sheet	7,680 Tips (80 Racks × 96 Tips)	S1120-1840-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de