



## Microdilution Tube System



### Product Highlights

- Microdilution Tube System offers a wide choice of configurations
- Ideal for automation / HTS and long term storage
- Reusable 96-place SBS racks containing individual, 8-strip or 12-strip tubes
- Temperature use -80 °C to 121 °C
- Tubes and racks are autoclavable
- Certified RNase, DNase, DNA and Pyrogen free
- Tubes and caps available separately.

## Product Applications

Multichannel and automated pipetting

Preparation of serial dilutions

Mixing

Storage and harvesting of cells

Cell growth assays

DNA screening

Long-term storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

|                       |                |
|-----------------------|----------------|
| Art. No.              | See variations |
| Autoclavable          | Yes            |
| Volume                | 1.2 ml         |
| Material              | Polypropylene  |
| DNA free              | Yes            |
| DNase free            | Yes            |
| Pyrogen free          | Yes            |
| RNase free            | Yes            |
| Operating temperature | -80°C to 121°C |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)















## More informations about Microdilution Tube System

### Organised storage!

Microdilution tube system comprising racks containing either individual, 8-strip or 12-strip 1.2 ml tubes. Ideal for serial dilutions, mixing, storage and harvesting of cells, cell growth assays, DNA screening and long-term storage.



## All Variations

| PRODUCT NAME   | PACKAGING SIZE            | ART. NO.   |
|--|---------------------------|------------|
| <br><a href="#">1.2 ml Reaction Tubes in Rack (96 individual Tubes)</a><br><b>Tubes per strip:</b> individual<br><b>Sterile:</b> No<br><br> <a href="#">Create PDF Data Sheet</a>            | 10 Pcs. (1 Box × 10 Pcs.) | E1750-0400 |
| <br><a href="#">1.2 ml Reaction Tubes in Rack (96 individual Tubes) (Sterile)</a><br><b>Tubes per strip:</b> individual<br><b>Sterile:</b> Yes<br><br> <a href="#">Create PDF Data Sheet</a> | 10 Pcs. (1 Box × 10 Pcs.) | E1750-0410 |
| <br><a href="#">1.2 ml Reaction Tubes in Rack (12 x 8-Strip Tubes)</a><br><b>Tubes per strip:</b> 8<br><b>Sterile:</b> No<br><br> <a href="#">Create PDF Data Sheet</a>                    | 10 Pcs. (1 Box × 10 Pcs.) | E1760-8400 |
| <br><a href="#">1.2 ml Reaction Tubes in Rack (12 x 8-Strip Tubes) (Sterile)</a><br><b>Tubes per strip:</b> 8<br><b>Sterile:</b> Yes<br><br> <a href="#">Create PDF Data Sheet</a>       | 10 Pcs. (1 Box × 10 Pcs.) | E1760-8410 |
| <br><a href="#">1.2 ml Reaction Tubes in Rack (8 x 12-Strip Tubes)</a><br><b>Tubes per strip:</b> 12<br><b>Sterile:</b> No<br><br> <a href="#">Create PDF Data Sheet</a>                 | 10 Pcs. (1 Box × 10 Pcs.) | E1770-9400 |
| <br><a href="#">1.2 ml Reaction Tubes in Rack (8 x 12-Strip Tubes) (Sterile)</a><br><b>Tubes per strip:</b> 12<br><b>Sterile:</b> Yes<br><br> <a href="#">Create PDF Data Sheet</a>      | 10 Pcs. (1 Box × 10 Pcs.) | E1770-9410 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)