



PlateOne® Deepwell Plates



Product Highlights

- Virgin polypropylene plates with SBS footprint
- Highly visible printed matrix
- The chemically resistant mats can be used with DMSO and other solvents
- Plate temperature use: -80 °C to 121 °C.
- Autoclavable
- Certified RNase, DNase, DNA and Pyrogen free
- Sterile plates are also pyrogen free

Product Applications

Sample preparation and storage
DNA and RNA isolation with high throughput
Plasmid purification with high throughput
Storage of cDNA or genomic banks
Storage of siRNA or oligonucleotide libraries
Protein analysis
Bacteria and yeast cultivation
Dilutions series

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

| | |
|-----------------------|----------------|
| Art. No. | See variations |
| Autoclavable | Yes |
| Material | Polypropylene |
| DNA free | Yes |
| DNase free | Yes |
| Human DNA free | Yes |
| RNase free | Yes |
| Color | Natural |
| Max. temp. | 121 °C |
| Temperature range | -80 to 121 °C |
| Temperature min. (SL) | -80 °C |
| Cut corner | A1 |
| g-Force max. | 5.000 x g |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about PlateOne® Deepwell Plates











Starlab 96- and 384-Deepwell Plates are high-performance plates for all manual and automated applications – from sample storage at -80 °C to DNA denaturation at 100 °C. Manufactured from pure polypropylene (PP), these plates offer innovative features: More dependability. More efficiency. More stability.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de















All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|---|---|------------|
|  0.5 ml PlateOne® 96 Deepwell Plate, Round Wells with Round Bottoms (Bulk) Volume: 0.5 ml Sterile: No Pack Size: 32 Plates (1 Box × 32 Plates)  Create PDF Data Sheet | 32 Plates (1 Box × 32 Plates) | S1896-5000 |
|  0.5 ml PlateOne® 96 Deepwell Plate, Round Wells with Round Bottoms (Sterile), Individually Wrapped Volume: 0.5 ml Sterile: Yes Pack Size: 32 Individually wrapped sachets (1 Box × 32 Individually wrapped sachets)  Create PDF Data Sheet | 32 Individually wrapped sachets (1 Box × 32 Individually wrapped sachets) | S1896-5110 |
|  0.5 ml PlateOne® 96 Deepwell Plate, Round Wells with Round Bottoms Volume: 0.5 ml Sterile: No Pack Size: 32 Plates (4 Packs × 8 Plates)  Create PDF Data Sheet | 32 Plates (4 Packs × 8 Plates) | S1896-5800 |
|  1.0 ml PlateOne® 96 Deepwell Plate, Round Wells with Round Bottoms (Bulk) Volume: 1.0 ml Sterile: No Pack Size: 20 Plates (1 Box × 20 Plates)  Create PDF Data Sheet | 20 Plates (1 Box × 20 Plates) | S1896-1000 |
|  1.0 ml PlateOne® 96 Deepwell Plate, Round Wells with Round Bottoms (Sterile), Individually Wrapped Volume: 1.0 ml Sterile: Yes Pack Size: 20 Individually wrapped sachets (1 Box × 20 Individually wrapped sachets)  Create PDF Data Sheet | 20 Individually wrapped sachets (1 Box × 20 Individually wrapped sachets) | S1896-1110 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de





| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|---|------------|
|  1.0 ml PlateOne® 96 Deepwell Plate, Round Wells with Round Bottoms Volume: 1.0 ml Sterile: No Pack Size: 20 Plates (4 Packs × 5 Plates)  Create PDF Data Sheet | 20 Plates (4 Packs × 5 Plates) | S1896-1800 |
|  2.0 ml PlateOne® 96 Deepwell Plate, Square Wells with V-Shape Bottoms (Bulk) Volume: 2.0 ml Sterile: No Pack Size: 20 Plates (1 Box × 20 Plates)  Create PDF Data Sheet | 20 Plates (1 Box × 20 Plates) | S1896-2000 |
|  2.0 ml PlateOne® 96 Deepwell Plate, Square Wells with V-Shape Bottoms (Sterile), Individually Wrapped Volume: 2.0 ml Sterile: Yes Pack Size: 20 Individually wrapped sachets (1 Box × 20 Individually wrapped sachets)  Create PDF Data Sheet | 20 Individually wrapped sachets (1 Box × 20 Individually wrapped sachets) | S1896-2110 |
|  2.0 ml PlateOne® 96 Deepwell Plate, Square Wells with V-Shape Bottoms Volume: 2.0 ml Sterile: No Pack Size: 20 Plates (4 Packs × 5 Plates)  Create PDF Data Sheet | 20 Plates (4 Packs × 5 Plates) | S1896-2800 |
|  200 µl PlateOne® 384 Deepwell Plate, Square Wells with V-Shape Bottoms (Bulk) Volume: 0.2 ml Sterile: No Pack Size: 40 Plates (1 Box × 40 Plates)  Create PDF Data Sheet | 40 Plates (1 Box × 40 Plates) | S1884-2000 |
|  200 µl PlateOne® 384 Deepwell Plate, Square Wells with V-Shape Bottoms Volume: 0.2 ml Sterile: No Pack Size: 40 Plates (4 Packs × 10 Plates)  Create PDF Data Sheet | 40 Plates (4 Packs × 10 Plates) | S1884-2400 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de






| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|--|------------|
|  <p>200 µl PlateOne® 384 Deepwell Plate, Square Wells with V-Shape Bottoms (Sterile), Individually Wrapped Volume: 0.2 ml Sterile: Yes Pack Size: 40 Individually wrapped sachets (1 Box × 40 Individually wrapped sachets)</p> <p> Create PDF Data Sheet</p> | 40 Individually wrapped sachets (1 Box × 40 Individually wrapped sachets) | S1884-2410 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



Accessories

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|---------------------------|------------|
|  Sealing Mat for 96-Well Plates with Square Wells, DMSO resistant | 24 Pcs. (1 Box × 24 Pcs.) | B1494-1124 |
|  Sealing Mat for 96-Well Plates with Round Wells, DMSO resistant | 24 Pcs. (1 Box × 24 Pcs.) | B1494-1224 |
|  Sealing Mat for 96-Well Plates with Square Wells, Piercable | 10 Pcs. (1 Box × 10 Pcs.) | E2896-1801 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de