



7.0 ml Bijou Containers, Sterile



Product Highlights

- Ultra clear sterile polystyrene container
- Centrifugation rate: 7.200 x g
- Temperature use: -10 °C to 70 °C
- 50 mm high x 23 mm diam (incl. cap)

Product Applications

Centrifugation

Sample preparation

Preparation of master mixes and buffers

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	See variations
Pack Size	700 Pcs. (1 Pack × 700 Pcs.)
Overall Dimensions	50 mm high x 23 mm diameter
Sterile	Yes
Autoclavable	No
Volume	7.0 ml
Material	Polystyrene
Max. temp.	70 °C
Temperature range	-10 to 70 °C
Min. temp.	-10 °C
g-Force max.	7.200 × g

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 7.0 ml Bijou Containers, Sterile





Ultra clear, sterile polystyrene bijou container with white polypropylene cap. Centrifugation rate 7200 x g. Temperature use: -10 °C to 70 °C.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	7.0 ml Polystyrene Bijou Container (Sterile) with Label Label Type: Plain Label	700 Pcs. (1 Pack × 700 Pcs.)	E1412-0711
 Create PDF Data Sheet			
	7.0 ml Polystyrene Bijou Container (Sterile) without Label Label Type: No Label	700 Pcs. (1 Pack × 700 Pcs.)	E1412-0710
 Create PDF Data Sheet			

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de