



Heat sealing films & foils



Product Highlights

- All seals are suitable for short-term storage and incubation.

Product Applications

Plate sealing

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

| | |
|----------|----------------|
| Art. No. | See variations |
| Seal | Heat sealing |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Heat sealing films & foils









All [heat seals](#) offered by STARLAB are suitable for PCR and one has superior optical clarity and is also suitable for Real-Time PCR.

Important note for Real-Time PCR (qPCR)

Please check that both the plate AND the sealing option are suitable for Real-Time PCR applications.



All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|---|-----------------------------|------------|
|  Peelable Film (PCR) Color: White Max. temp.: 90 °C  Create PDF Data Sheet | 100 Pcs. (1 Box × 100 Pcs.) | E2796-3060 |
|  Thermo Seal Film (PCR) Color: Silver Max. temp.: 120 °C  Create PDF Data Sheet | 100 Pcs. (1 Box × 100 Pcs.) | E2796-3070 |
|  Easy Pierce Film (PCR) Color: Silver Max. temp.: 80 °C  Create PDF Data Sheet | 100 Pcs. (1 Box × 100 Pcs.) | E2796-3090 |
|  Optically Clear Film (qPCR) Color: Clear Max. temp.: 120 °C  Create PDF Data Sheet | 100 Pcs. (1 Box × 100 Pcs.) | E2796-3020 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de