



Cell Strainers (Sterile)



Product Highlights

- Uniform nylon mesh size of 40 μm , 70 μm , and 100 μm (colour-coded)
- Extended lip on the strainer enables aseptic handling with forceps
- Designed to fit into a 50 ml conical centrifuge tube
- RNase, DNase and Pyrogen free
- Gamma sterilised (SAL 10^{-6})

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	See variations
Pack Size	50 Pcs. (50 Box × 1 Pcs.)
Sterile	Yes
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes
Non-cytotoxic	Yes
For use with (SL)	50 ml centrifuge tube

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Cell Strainers (Sterile)







STARLAB Cell Strainers are intended as a fast, easy and simple alternative to traditional filtration when dissociating stem cells and other tissue-derived primary cells. They yield consistently uniform single-cell suspensions that are ideal for preparation of samples for flow cytometry, Fluorescence-Activated Cell Sorting (FACS), and other applications related to cell separation. Available in 40 µm, 70 µm, and 100 µm color-coded mesh variants. Cell strainers are made from USP Class VI materials.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	40 µm Cell Strainer, Blue (Sterile) Pore size: 40 µm	50 Pcs. (50 Box × 1 Pcs.)	CC8111-0042
	 Create PDF Data Sheet		
	70 µm Cell Strainer, White (Sterile) Pore size: 70 µm	50 Pcs. (50 Box × 1 Pcs.)	CC8111-0072
	 Create PDF Data Sheet		
	100 µm Cell Strainer, Yellow (Sterile) Pore size: 100 µm	50 Pcs. (50 Box × 1 Pcs.)	CC8111-0102
	 Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de