



15 ml Centrifuge Tubes



Product Highlights

- Clear, sterile polypropylene tubes with green HDPE screw caps
- White writing surface and clear, printed graduations
- Perfect for centrifugation of cells: USP Class VI, non-cytotoxic and sterile (SAL 10^{-6}) material
- Certified RNase, DNase, DNA, heavy-metal and Pyrogen free

Product Applications

Centrifugation

Sample preparation

Preparation of master mixes and buffers

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	See variations
Dimensions (W × H)	22 mm × 120.5 mm
External width	22 mm
Sterile	Yes
Autoclavable	Yes
Volume	15 ml
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Free of heavy metals	Yes
Human DNA free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Clear
g-Force max.	12,000 x g
Non-cytotoxic	Yes
Tube Base	Conical

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de







More informations about 15 ml Centrifuge Tubes

STARLAB's new centrifuge tubes are designed for the centrifugation of cells, but can also be used for standard centrifugation, preparation and storage applications. USP Class VI polymers and certified non-toxicity ensure the viability and integrity of cells. The sterility of assurance level (SAL) of 10^{-6} fulfills highest demands on sterility, meeting EN 556-1 standard (Sterilisation of medical devices). Racked tubes are delivered in recyclable cardboard racks.

Maximum centrifugation rates: Conical Tubes 15 & 50 ml: 12,000 x g and Skirted Tube 50 ml: 6,000 x g



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	15 ml Centrifuge Tube, Conical, Loose (Sterile) Packaging Unit: Bags	500 Tubes (20 Bags × 25 Tubes)	E1415-0100
 Create PDF Data Sheet			
	15 ml Centrifuge Tube, Conical, Racked (Sterile) Packaging Unit: Racks	500 Tubes (20 Racks × 25 Tubes)	E1415-0700
 Create PDF Data Sheet			

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de