




Flat Gel-Loading Tips



Product Highlights

- Flexible, long and narrow
- Loading samples onto gels (eg. polyacrylamide or agarose)
- Handling of small volumes
- Compatible with a wide range of pipettes
- Different pack types available
- Certified

Product Applications

Gel electrophoresis

Loading samples onto gels (eg. polyacrylamide or agarose)

Handling of small volumes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	See variations
Sterile	No
Volume	10 µl
Pack Size	200 Tips (1 Rack × 200 Tips)
Color	Natural
Tip Type	flat gel-loading

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Flat Gel-Loading Tips

All tips are certified RNase and DNase free. Packaging of these tips differs from the TipOne® range.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 10 µl Flat Gel-Loading Tip (45 mm long and 0.20 mm thick) Tip Thickness: 0.20 mm  Create PDF Data Sheet	200 Tips (1 Rack × 200 Tips)	I1012-2600
 10 µl Flat Gel-Loading Tip (45 mm long and 0.40 mm thick) Tip Thickness: 0.40 mm  Create PDF Data Sheet	200 Tips (1 Rack × 200 Tips)	I1012-4600

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de