



## 200 µl Bevelled TipOne® Tip, Natural



### Product Highlights

- Universal fit. One tip fits all! The tip collar is engineered to fit all popular pipettes.
- Highly polished moulds for an extremely smooth surface to reduce sample retention
- Certified RNase, DNase, DNA and Pyrogen free.

## Product Applications

Pipetting liquids

Dispensing liquids

Mixing liquids

Filling plates and reaction vessels

Phase extraction and removal of supernatants

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	See variations
Autoclavable	Yes
Filter tip	Non-filter
Volume range	1.0 - 200 µl
Volume	200 µl
Tip Type	bevelled
Material	Polypropylene
Graduations	Yes
Color	Natural

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about 200 µl Bevelled TipOne® Tip, Natural

### Improved by you - New System Components Overview

#### The perfect rack

- › Simple! The system comprises only two racks (one small, one large) for less complexity
- › Small rack for use with tip volumes 10 µl to 300 µl; large rack for use with tip volumes 1000 µl and 1250µl XL
- › Quick and easy identification of tips by colour-coded wafers with laser printed information and transparent lid
- › Stackable
- › Reusable
- › Autoclavable

#### The clever refill

- › The new refill is divided into two separate, individually-wrapped stacks of five wafers
- › Completely sealed, the tips are safe from contamination by the tight fitting top shell and unique spacers
- › Economical, space-saving design
- › Autoclavable when racked

#### The self-standing bag

- › New, self-standing bag for easy and convenient access and storage
- › Tamper-proof, easy-open perforation strip
- › Resealable zip-lock for contamination-free storage
- › Recyclable polyethylene bag
- › Economical – ideal for teaching labs or non-critical work

#### The compact stack rack

- › Next generation design with wide base for great stability
- › Transparent lid covers the tips when not in use
- › 960 tips for high throughput use, ideal for multichannel users
- › The complete stack rack is autoclavable and recyclable

#### The unique filter refill













- › A brand new, economical and ecofriendly TipOne® filter tip refill system
- › Fast, easy and safe: The filter tips are refilled with the top and bottom shell intact, thus providing a fast and completely protected refill system that will not compromise sterility
- › Using sterile filter refills to reload racks directly during periods of long use in hoods and PCR-workstations saves valuable space
- › Not only does the refill system generate less waste, the bottom shells are nestable so your waste takes up less storage space too

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)















## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">200 µl Bevelled TipOne® Tip, Natural, Racks (non-sterile)</a> <b>Sterile:</b> No <b>Pack Size:</b> 960 Tips (10 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Racks × 96 Tips)	S1111-1800
 <a href="#">200 µl Bevelled TipOne® Tip, Natural, Racks (non-sterile), Case</a> <b>Sterile:</b> No <b>Pack Size:</b> 7,680 Tips (80 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	7,680 Tips (80 Racks × 96 Tips)	S1111-1800-C
 <a href="#">200 µl Bevelled TipOne® Tip, Natural, Racks (sterile)</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 960 Tips (10 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Racks × 96 Tips)	S1111-1810
 <a href="#">200 µl Bevelled TipOne® Tip, Natural, Racks (sterile), Case</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 7,680 Tips (80 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	7,680 Tips (80 Racks × 96 Tips)	S1111-1810-C
 <a href="#">200 µl Bevelled TipOne® Tip, Natural, Refills (non-sterile)</a> <b>Sterile:</b> No <b>Pack Size:</b> 960 Tips (10 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Refills × 96 Tips)	S1111-1700
 <a href="#">200 µl Bevelled TipOne® Tip, Natural, Refills (non-sterile), Case</a> <b>Sterile:</b> No <b>Pack Size:</b> 9,600 Tips (100 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	9,600 Tips (100 Refills × 96 Tips)	S1111-1700-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)







PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">200 µl Bevelled TipOne® Tip, Natural, Stack rack (non-sterile)</a> <b>Sterile:</b> No <b>Pack Size:</b> 960 Tips (1 Stack rack × 960 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (1 Stack rack × 960 Tips)	S1111-1200
 <a href="#">200 µl Bevelled TipOne® Tip, Natural, Stack racks (non-sterile), Case</a> <b>Sterile:</b> No <b>Pack Size:</b> 9,600 Tips (10 Stack racks × 960 Tips)  <a href="#">Create PDF Data Sheet</a>	9,600 Tips (10 Stack racks × 960 Tips)	S1111-1200-C
 <a href="#">200 µl Bevelled TipOne® Tip, Natural, Stack rack (sterile)</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 960 Tips (1 Stack rack × 960 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (1 Stack rack × 960 Tips)	S1111-1210
 <a href="#">200 µl Bevelled TipOne® Tip, Natural, Stack racks (sterile), Case</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 9,600 Tips (10 Stack racks × 960 Tips)  <a href="#">Create PDF Data Sheet</a>	9,600 Tips (10 Stack racks × 960 Tips)	S1111-1210-C
 <a href="#">200 µl Bevelled TipOne® Tip, Natural, Bag (non-sterile)</a> <b>Sterile:</b> No <b>Pack Size:</b> 1,000 Tips (1 Bag × 1000 Tips)  <a href="#">Create PDF Data Sheet</a>	1,000 Tips (1 Bag × 1000 Tips)	S1111-1000
 <a href="#">200 µl Bevelled TipOne® Tip, Natural, Bags (non-sterile), Case</a> <b>Sterile:</b> No <b>Pack Size:</b> 25,000 Tips (25 Bags × 1000 Tips)  <a href="#">Create PDF Data Sheet</a>	25,000 Tips (25 Bags × 1000 Tips)	S1111-1000-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



PRODUCT NAME		PACKAGING SIZE	ART. NO.
	<a href="#">200 µl Bevelled TipOne® Tip, Natural, Refills (sterile)</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 960 Tips (10 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Refills × 96 Tips)	S1111-1710
	<a href="#">200 µl Bevelled TipOne® Tip, Natural, Refills (sterile), Case</a> <b>Sterile:</b> Yes <b>Pack Size:</b> 7,680 Tips (80 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	7,680 Tips (80 Refills × 96 Tips)	S1111-1710-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)