



## Ultra Clear, High Speed Microcentrifuge Tubes



### Product Highlights

- Ultra clear polypropylene tubes
- Firm and positive cap closure
- Suitable for boiling
- Rated to over 20,000 xg
- Certified RNase, DNase, DNA and PCR Inhibitor free

## Product Applications

Sample preparation

Centrifugation

Sample storage

Sample boiling or freezing

Certified RNase, DNase, DNA and PCR Inhibitor free

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	See variations
Autoclavable	Yes
Volume	1.5 ml
Pack Size	500 Tubes (1 Bag × 500 Tubes)
Material	Polypropylene
DNA free	Yes
DNase free	Yes
PCR inhibitor free	Yes
RNase free	Yes
Color	Various

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about Ultra Clear, High Speed Microcentrifuge Tubes









Polypropylene tubes with smaller lid footprint for better fit inside [sample storage](#) boxes. Amber tubes for light sensitive samples. Rated to over 20,000 x g. Temp. range: -80 °C to 121 °C. Autoclavable. Certified RNase, DNase, DNA and PCR Inhibitor free.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">1.5 ml Homopolymer Microcentrifuge Tube, Amber</a> <b>Color:</b> Amber  <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	I1415-1007
 <a href="#">1.5 ml Homopolymer Microcentrifuge Tube, Natural</a> <b>Color:</b> Natural  <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	I1415-1000
 <a href="#">1.5 ml Homopolymer Microcentrifuge Tube, Mixed</a> <b>Color:</b> Mixed  <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	I1415-1008
 <a href="#">1.5 ml Homopolymer Microcentrifuge Tube, Violet</a> <b>Color:</b> Violet  <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	I1415-1005

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## Accessories

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">Cap-Locks for 1.5 / 2.0 ml Microcentrifuge Tubes</a>	100 Pcs. (1 Bag × 100 Pcs.)	I1415-1508

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)