



Bottle Top Filtration Unit, Full Assembly, 0.22 µm, 500 ml



Product Highlights

- 0.22 µm PES membrane
- High speed filtration
- Funnels compatible with all common bottles (GL 45 compatible)
- Temperature range: -20°C to +50 C°
- Sterile SAL10⁻⁶ (ISO11137)
- Certified RNase, DNase free and non-pyrogenic
- Wide bottle base for stability
- Wide funnel top for fast and efficient use
- Hose connector designed to fit multiple hose diameters

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	CC8221-5226
Sterile	Yes
Volume	500 ml
Pack Size	12 Pcs. (1 Box × 12 Pcs.)
Free of endotoxins	Yes
Free of heavy metals	Yes
Pyrogen free	Yes
RNase free	Yes
Non-cytotoxic	Yes
Pore size	0.22 µm
System	Full Assembly

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Bottle Top Filtration Unit, Full Assembly, 0.22 µm, 500 ml

OPEN WIDE!













Bottle Top Filtration units with a wide base for stability and an easily accessible, wide funnel top for efficient use. The clever hose connector is designed to fit multiple hose diameters. The media bottle caps have large knurls for easy loosening or tightening. Each product is manufactured in a Class 100,000 cleanroom environment and is printed with the batch number and expiry date for easy traceability.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de







All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Bottle Top Filtration Unit, Full Assembly, 0.22 µm, 150 ml Pore size: 0.22 µm Volume: 150 ml  Create PDF Data Sheet	12 Pcs. (1 Box × 12 Pcs.)	CC8221-3226
 Bottle Top Filtration Unit, Full Assembly, 0.22 µm, 250 ml Pore size: 0.22 µm Volume: 250 ml  Create PDF Data Sheet	12 Pcs. (1 Box × 12 Pcs.)	CC8221-4226
 Bottle Top Filtration Unit, Full Assembly, 0.22 µm, 500 ml Pore size: 0.22 µm Volume: 500 ml  Create PDF Data Sheet	12 Pcs. (1 Box × 12 Pcs.)	CC8221-5226
 Bottle Top Filtration Unit, Full Assembly, 0.22 µm, 1000 ml Pore size: 0.22 µm Volume: 1000 ml  Create PDF Data Sheet	12 Pcs. (1 Box × 12 Pcs.)	CC8221-6226
 Bottle Top Filtration Unit, Full Assembly, 0.45 µm, 150 ml Pore size: 0.45 µm Volume: 150 ml  Create PDF Data Sheet	12 Pcs. (1 Box × 12 Pcs.)	CC8221-3456
 Bottle Top Filtration Unit, Full Assembly, 0.45 µm, 250 ml Pore size: 0.45 µm Volume: 250 ml  Create PDF Data Sheet	12 Pcs. (1 Box × 12 Pcs.)	CC8221-4456

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Bottle Top Filtration Unit, Full Assembly, 0.45 µm, 500 ml Pore size: 0.45 µm Volume: 500 ml  Create PDF Data Sheet	12 Pcs. (1 Box × 12 Pcs.)	CC8221-5456
 Bottle Top Filtration Unit, Full Assembly, 0.45 µm, 1000 ml Pore size: 0.45 µm Volume: 1000 ml  Create PDF Data Sheet	12 Pcs. (1 Box × 12 Pcs.)	CC8221-6456

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de