



## Adapter for pipettes



### Product Highlights

- Spare part for use with Aspirator's Hand Controller
- For use with Pasteur, aspiration and graduated pipettes
- For pipettes with an outer diameter of 5.5 - 8.0 mm
- Autoclavable

## Product Applications

Aspiration of Liquids

Cell-Culture

Sterile Filtration

De-gassing of buffers

Solidphase Extraction

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	N2400-9016
Pack Size	1 Piece (1 Box × 1 Piece)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about Adapter for pipettes

This is a spare part for the [Aspirator's](#) Hand Controller (N2400-9001). Adapter fits Pasteur, aspiration and graduated [pipettes](#) with an outer diameter of 5.5 - 8.0 mm. An illustrated description of how to replace the adapter can be found in the manual for the hand controller (see download area).