



## Plate Sealing



### Product Highlights

- Provides a secure seal
- Use for short-term storage, incubation and containment of biohazards
- 50 µm clear polyester seal with acrylic adhesive
- White poly-coated paper release liner
- Suitable for ELISA, EIA and similar assays
- Eliminates the edge effect in sensitive ELISA assays
- Temperature range: -40°C to 120°C
- Dimensions: 146 x 79.4 mm

## Product Applications

Plate Sealing

Suitable for ELISA, EIA and similar assays

Eliminates the edge effect in sensitive ELISA assays

Use for short-term storage, incubation and containment of biohazards

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	See variations
Pack Size	100 Pcs. (1 Box × 100 Pcs.)
Material	Polyester
DNase free	No
Free of endotoxins	No
Pyrogen free	No
RNase free	No
Color	Clear
Seal	Self-adhesive

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)







## More informations about Plate Sealing

A sterile and non-sterile variant of self-adhesive, clear polyester film which is suitable for ELISA, EIA and similar assays. Manufactured from 50 µm polyester, these non-piercable seals provide a secure seal for every well to eliminate edge effects in sensitive ELISA assays. Also suitable for use with Tissue Culture plates for short-term storage, incubation, and containment of biohazards.



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">Polyester Sealing Film, Clear</a> <b>Sterile:</b> No  <a href="#">Create PDF Data Sheet</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-0101
 <a href="#">Polyester Sealing Film, Clear</a> <b>Sterile:</b> Yes  <a href="#">Create PDF Data Sheet</a>	100 Pcs. (1 Box × 100 Pcs.)	E2796-0102

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)