



## Tips for Tecan® & Packard® - Conductive (Non-Filter)



### Product Highlights

- All robotic tips are certified RNase, DNase and DNA free
- No special set-up required!
- Tips for Tecan® & Packard®

## Product Applications

PCR Setup

Assay-Setup


Serial Dilutions

Distributing Reagents

Sample or Reagent Transfer

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	See variations
Pack Size	2,304 Tips (24 Racks × 96 Tips)
Sterile	No
Filter tip	Non-filter
Rack Style	White Hanging Rack
Conductive	Yes
Head	Varispan with 4 or 8 head probe
For use with	Tecan® & Packard®

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about Tips for Tecan® & Packard® - Conductive (Non-Filter)

STARLAB robotic tips are engineered to the same specifications as the original equipment manufacturers' tips to ensure compatibility and a consistent fit. Don't risk having your time and results compromised with defective materials and poor performance! With STARLAB, you can put your trust in disposables from a company that offers much more than just tips.

All robotic tips are certified RNase, DNase and DNA free

### Quality control

Quality control is the essential factor which gives you the freedom to walk away from your automated workstation and get on with other tasks. You need the assurance that the tip will not let you down. STARLAB robotic tips are put through stringent and thorough testing. Every lot of tips is subject to an array of demanding tests, from the geometry of the tip itself to automation-specific requirements such as clean ejection from the liquid handling head. Every lot is required to meet or exceed the parameters necessary to give you that confidence to walk away.

Tips are tested on the original robotic equipment so there is no doubt they will perform to the same requirements as the disposables from the robot manufacturer themselves.

### No special set-up required!





STARLAB robotic tips require no instrument hardware or software adjustment, thereby providing a simple, safe and direct high-performance substitute.

### Purchasing flexibility

STARLAB likes to make things as simple as possible for its customers. We realise that you may need to order robotic consumables in large quantities, but may not have the space to store them. STARLAB can help with call-off orders where you buy and own the product, but we store it in our warehouse and 'call off' delivery to you, as and when you need it. Alternatively, you may place a standing order and have a predetermined quantity delivered at regular intervals. Contact us to discuss your options.



## All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	<a href="#">50µl Conductive Tip for Tecan® + Packard® (non-sterile)</a> <b>Volume:</b> 50 µl	2,304 Tips (24 Racks × 96 Tips)	R1181-5709
	 <a href="#">Create PDF Data Sheet</a>		
	<a href="#">200µl Conductive Tip for Tecan® &amp; Packard® (non-sterile)</a> <b>Volume:</b> 200 µl	2,304 Tips (24 Racks × 96 Tips)	R1181-0709
	 <a href="#">Create PDF Data Sheet</a>		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)