



## TwisStar Rack



### Product Highlights

- Polypropylene rack
- Rotate and lock together to achieve various configurations

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	E2344-1099
Dimensions (W × D × H)	257 mm × 108 mm × 108 mm
Dimensions (W × D)	257 mm × 108 mm
Autoclavable	Yes
Pack Size	1 Piece (1 Rack × 1 Piece)
Rack Type	Multi-Tube Rack
Material	Polypropylene
Color	Mixed
Max. centrifuge tubes or screwcap microtubes	50 ml tubes (max. 30mm $\varnothing$ ) x 2 per section 15 ml tubes (max. 18mm $\varnothing$ ) x 8 per section 1.5 / 2.0 ml tubes (max. 11mm $\varnothing$ ) x 12 per section 0.5 ml tubes (max. 8.0 mm $\varnothing$ ) x 6 per section PCR Plate 1 plate across 2 sections 8-Tube PCR Stri

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about TwisStar Rack

TwisStar – polypropylene rack comprising four sections that can be rotated and locked together to achieve various configurations for holding different tubes and a PCR plate.

Dimensions (WxDxH): 257 x 108 x 108 mm (combined) and 69 x 108 x 108 mm (individual section)