



Media Bottle, PS, 250 ml (sterile)



Product Highlights

- 250 ml PS Media Bottle
- Compatible with all common Bottle Top Filters (GL 45 compatible)
- Wide base for stability
- Temperature range: -20 °C to +50 °C
- Sterile SAL10⁻⁶ (ISO1137)
- Certified RNase, DNase free and non-pyrogenic
- Wide bottle top for fast and efficient use
- Ergonomically designed for easier handling to help prevent spills.

Product Applications

Sterile filtration
Media preparation
Media storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

| | |
|----------------------|---------------------------|
| Art. No. | CC8221-4000 |
| Sterile | Yes |
| Autoclavable | No |
| Volume | 250 ml |
| Pack Size | 24 Pcs. (1 Box × 24 Pcs.) |
| Material | Polystyrene |
| DNase free | Yes |
| Free of endotoxins | Yes |
| Free of heavy metals | Yes |
| Pyrogen free | Yes |
| RNase free | Yes |
| Non-cytotoxic | Yes |
| System | Media Bottle |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Media Bottle, PS, 250 ml (sterile)

OPEN WIDE!









Media Bottles with a wide base for stability and a wide, easily accessible bottle top for efficient use. The reservoir bottle caps have large knurls for easy loosening or tightening. Each product is manufactured in a Class 100,000 cleanroom environment and is printed with the batch number and expiry date for easy traceability.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|---|---------------------------|-------------|
|  Media Bottle, PS, 150 ml (sterile) Volume: 150 ml  Create PDF Data Sheet | 24 Pcs. (1 Box × 24 Pcs.) | CC8221-3000 |
|  Media Bottle, PS, 250 ml (sterile) Volume: 250 ml  Create PDF Data Sheet | 24 Pcs. (1 Box × 24 Pcs.) | CC8221-4000 |
|  Media Bottle, PS, 500 ml (sterile) Volume: 500 ml  Create PDF Data Sheet | 24 Pcs. (1 Box × 24 Pcs.) | CC8221-5000 |
|  Media Bottle, PS, 1000 ml (sterile) Volume: 1000 ml  Create PDF Data Sheet | 24 Pcs. (1 Box × 24 Pcs.) | CC8221-6000 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de