



## Bottle Top Filtration Unit, Funnel Only, 0.45 µm, 250 ml (sterile)



### Product Highlights

- 0.45 µm PES membrane
- High speed filtration
- Compatible with all common bottles (GL 45 compatible)
- Temperature range: -20 °C to +50 °C
- Sterile SAL10<sup>-6</sup> (ISO1137) / 3-year shelf life
- Certified RNase, DNase free and non-pyrogenic
- Hose connector designed to fit multiple hose diameters

## Product Applications

Sterile filtration  
Media preparation

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	CC8222-4456
Pack Size	24 Pcs. (1 Box × 24 Pcs.)
Sterile	Yes
Autoclavable	No
Volume	250 ml
DNase free	Yes
Free of endotoxins	Yes
Free of heavy metals	Yes
Pyrogen free	Yes
RNase free	Yes
Non-cytotoxic	Yes
Pore size	0.45 µm
System	Funnel Only

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about Bottle Top Filtration Unit, Funnel Only, 0.45 µm, 250 ml (sterile)













Bottle Top Filtration Funnels with PES membrane for fast filtration. The clever hose connector is designed to fit multiple hose diameters. Each product is manufactured in a Class 100,000 cleanroom environment and is printed with the batch number and expiry date for easy traceability.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)







## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">Bottle Top Filtration Unit, Funnel Only, 0.22 µm, 150 ml (sterile)</a> <b>Pore size:</b> 0.22 µm <b>Volume:</b> 150 ml   <a href="#">Create PDF Data Sheet</a>	24 Pcs. (1 Box × 24 Pcs.)	CC8222-3226
 <a href="#">Bottle Top Filtration Unit, Funnel Only, 0.22 µm, 250 ml (sterile)</a> <b>Pore size:</b> 0.22 µm <b>Volume:</b> 250 ml   <a href="#">Create PDF Data Sheet</a>	24 Pcs. (1 Box × 24 Pcs.)	CC8222-4226
 <a href="#">Bottle Top Filtration Unit, Funnel Only, 0.22 µm, 500 ml (sterile)</a> <b>Pore size:</b> 0.22 µm <b>Volume:</b> 500 ml   <a href="#">Create PDF Data Sheet</a>	24 Pcs. (1 Box × 24 Pcs.)	CC8222-5226
 <a href="#">Bottle Top Filtration Unit, Funnel Only, 0.22 µm, 1000 ml (sterile)</a> <b>Pore size:</b> 0.22 µm <b>Volume:</b> 1000 ml   <a href="#">Create PDF Data Sheet</a>	24 Pcs. (1 Box × 24 Pcs.)	CC8222-6226
 <a href="#">Bottle Top Filtration Unit, Funnel Only, 0.45 µm, 150 ml (sterile)</a> <b>Pore size:</b> 0.45 µm <b>Volume:</b> 150 ml   <a href="#">Create PDF Data Sheet</a>	24 Pcs. (1 Box × 24 Pcs.)	CC8222-3456
 <a href="#">Bottle Top Filtration Unit, Funnel Only, 0.45 µm, 250 ml (sterile)</a> <b>Pore size:</b> 0.45 µm <b>Volume:</b> 250 ml   <a href="#">Create PDF Data Sheet</a>	24 Pcs. (1 Box × 24 Pcs.)	CC8222-4456

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



PRODUCT NAME		PACKAGING SIZE	ART. NO.
	<a href="#">Bottle Top Filtration Unit, Funnel Only, 0.45 µm, 500 ml (sterile)</a> <b>Pore size:</b> 0.45 µm <b>Volume:</b> 500 ml	24 Pcs. (1 Box × 24 Pcs.)	CC8222-5456
	 <a href="#">Create PDF Data Sheet</a>		
	<a href="#">Bottle Top Filtration Unit, Funnel Only, 0.45 µm, 1000 ml (sterile)</a> <b>Pore size:</b> 0.45 µm <b>Volume:</b> 1000 ml	24 Pcs. (1 Box × 24 Pcs.)	CC8222-6456
	 <a href="#">Create PDF Data Sheet</a>		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)