



Filtration Funnels (sterile)



Product Highlights

- Pore sizes 0.22 and 0.45 µm
- High speed filtration
- Compatible with all common bottles (GL 45 compatible)
- Temperature range: -20 °C to +50 °C
- Sterile SAL10⁻⁶ (ISO1137)
- Certified RNase, DNase free and non-pyrogenic
- Hose connector designed to fit multiple hose diameters

Product Applications

Sterile filtration
Media preparation

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

| | |
|----------------------|---------------------------|
| Art. No. | See variations |
| Sterile | Yes |
| Autoclavable | No |
| Pack Size | 24 Pcs. (1 Box × 24 Pcs.) |
| DNase free | Yes |
| Free of endotoxins | Yes |
| Free of heavy metals | Yes |
| Pyrogen free | Yes |
| RNase free | Yes |
| Non-cytotoxic | Yes |
| System | Funnel Only |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Filtration Funnels (sterile)













Bottle Top Filtration Funnels with PES mebrane for fast filtration. The clever hose connector is designed to fit multiple hose diameters. Each product is manufactured in a Class 100,000 cleanroom environment and is printed with the batch number and expiry date for easy traceability.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de







All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|---------------------------|-------------|
|  Bottle Top Filtration Unit, Funnel Only, 0.22 µm, 150 ml (sterile) Pore size: 0.22 µm Volume: 150 ml  Create PDF Data Sheet | 24 Pcs. (1 Box × 24 Pcs.) | CC8222-3226 |
|  Bottle Top Filtration Unit, Funnel Only, 0.22 µm, 250 ml (sterile) Pore size: 0.22 µm Volume: 250 ml  Create PDF Data Sheet | 24 Pcs. (1 Box × 24 Pcs.) | CC8222-4226 |
|  Bottle Top Filtration Unit, Funnel Only, 0.22 µm, 500 ml (sterile) Pore size: 0.22 µm Volume: 500 ml  Create PDF Data Sheet | 24 Pcs. (1 Box × 24 Pcs.) | CC8222-5226 |
|  Bottle Top Filtration Unit, Funnel Only, 0.22 µm, 1000 ml (sterile) Pore size: 0.22 µm Volume: 1000 ml  Create PDF Data Sheet | 24 Pcs. (1 Box × 24 Pcs.) | CC8222-6226 |
|  Bottle Top Filtration Unit, Funnel Only, 0.45 µm, 150 ml (sterile) Pore size: 0.45 µm Volume: 150 ml  Create PDF Data Sheet | 24 Pcs. (1 Box × 24 Pcs.) | CC8222-3456 |
|  Bottle Top Filtration Unit, Funnel Only, 0.45 µm, 250 ml (sterile) Pore size: 0.45 µm Volume: 250 ml  Create PDF Data Sheet | 24 Pcs. (1 Box × 24 Pcs.) | CC8222-4456 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|---------------------------|-------------|
|  Bottle Top Filtration Unit, Funnel Only, 0.45 µm, 500 ml (sterile) Pore size: 0.45 µm Volume: 500 ml  Create PDF Data Sheet | 24 Pcs. (1 Box × 24 Pcs.) | CC8222-5456 |
|  Bottle Top Filtration Unit, Funnel Only, 0.45 µm, 1000 ml (sterile) Pore size: 0.45 µm Volume: 1000 ml  Create PDF Data Sheet | 24 Pcs. (1 Box × 24 Pcs.) | CC8222-6456 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de