




Storage box 5.0ml tubes, blue



Product Highlights

- 25-place storage box for 5.0 ml Centrifuge Tubes
- Transparent lid with printed reference grid for easier sample identification
- One-way keyed lid prevents sample mix-up
- Air vents minimize condensation
- Base with drain holes
- Stackable
- Temperature use: -90°C to 80°C
- Dimensions (WxDxH): 133 x 133 x 67 mm (with lid)
- Cover made of polystyrene / Base made of high impact polystyrene

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

| | |
|------------------------|--------------------------|
| Art. No. | E2325-5021 |
| Dimensions (W × D × H) | 133 mm × 133 mm × 67 mm |
| Dimensions (W × D) | 133 mm × 133 mm |
| Autoclavable | No |
| Pack Size | 3 Pcs. (1 Pack × 3 Pcs.) |
| Material | Polystyrene |
| Color | Blue |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Storage box 5.0ml tubes, blue

25-place stackable, storage box for your 5.0 ml Centrifuge Tubes. The transparent lid features a printed reference grid for easy sample identification. Air vents minimise condensation and drain holes in the base aid lid removal following freezing. Suitable for use down to -90°C .

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de