



## Rack 5.0ml tubes, blue



### Product Highlights

- 25-place rack for 5.0 ml Centrifuge Tubes
- Also holds 15 ml Centrifuge Tubes
- Numbered wells to aide identification
- Temperature use: -90°C to 121°C
- Autoclavable
- Dimensions (WxDxH): 146.8 x 113.8 x 38 mm
- Made of Polypropylene

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	E2325-5011
Pack Size	5 Pcs. (1 Pack × 5 Pcs.)
Dimensions (W × D × H)	146.8 mm × 113.8 mm × 38 mm
Dimensions (W × D)	146.8 mm × 113.8 mm
External width	146.8 mm
Rack Type	Tube Rack
Material	Polypropylene
Color	Blue

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about Rack 5.0ml tubes, blue

25-place polypropylene rack for 5.0 ml and 15 ml concial Centrifuge Tubes. The rack features numbered wells to aide sample identification and organisation on the bench. The rack is reusable and autoclavable, and can also be used to store samples in the freezer.