



## Centrifuge tube, 5ml, flip cap, conical, amber



### Product Highlights

- RNase, DNase, DNA & Pyrogen free
- Temperature use: -196°C to 121°C
- Autoclavable
- Centrifugation rate: 25,000 x g
- Molded graduations at 0.5 ml increments
- Etched writing surface on flat cap and side wall
- Cap center with "membrane" area for easy puncture with a needle
- Dimensions (HxØ): 60.2 x 15.8 mm
- Made of ultra-clear polypropylene

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	E1450-5007
Pack Size	200 Tubes (2 Bags × 100 Tubes)
Sterile	No
Volume	5.0 ml
Cap style	Flip Cap
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Amber
g-Force max.	25.000 x g

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about Centrifuge tube, 5ml, flip cap, conical, amber

Do you need to prepare samples up to 5.0 ml but don't want to use 15 ml tubes? Save on plastic waste, with our 5.0 ml Centrifuge Tubes!











- Ideal for biological or chemical solution preparation
- Easier sample access lowers the risk of contamination
- Saves on space and plastic, compared to 15 ml tubes
- Suitable for microbial incubation and pellet preparation
- Choice of cap styles available - flip top or screw cap
- Sterile options available
- Amber tube for light sensitive reactions
- Various accesories available

No processing additives such as slip agents, biocides, plasticizers or mold release agents are used during production.

Boil-proof design : When the internal pressure approaches an unsafe level, the cap partially opens, softly, in a controlled manner, releasing excess pressure, but preventing any splashing of contents. In applications over 80 °C, we recommend cap lock clips.



## All Variations






PRODUCT NAME		PACKAGING SIZE	ART. NO.
	<a href="#">Centrifuge tube, 5ml, flip cap, conical, natural</a> <b>Cap style:</b> Flip Cap <b>Sterile:</b> No <b>Color:</b> Natural   <a href="#">Create PDF Data Sheet</a>	200 Tubes (2 Bags × 100 Tubes)	E1450-5000
	<a href="#">Centrifuge tube, 5ml, flip cap, conical, amber</a> <b>Cap style:</b> Flip Cap <b>Sterile:</b> No <b>Color:</b> Amber   <a href="#">Create PDF Data Sheet</a>	200 Tubes (2 Bags × 100 Tubes)	E1450-5007
	<a href="#">Centrifuge tube, 5ml, flip cap, conical, natural, sterile</a> <b>Cap style:</b> Flip Cap <b>Sterile:</b> Yes <b>Color:</b> Natural   <a href="#">Create PDF Data Sheet</a>	200 Tubes (10 Bags × 20 Tubes)	E1450-5010
	<a href="#">Centrifuge tube, 5ml, screw tube only, conical, natural</a> <b>Cap style:</b> Screw Cap <b>Sterile:</b> No <b>Color:</b> Natural   <a href="#">Create PDF Data Sheet</a>	200 Tubes (2 Bags × 100 Tubes)	E1450-5040
	<a href="#">Centrifuge tube, 5.0ml, with standard screw cap, conical, natural white, sterile</a> <b>Cap style:</b> Screw Cap <b>Sterile:</b> Yes <b>Color:</b> Natural   <a href="#">Create PDF Data Sheet</a>	200 Tubes (10 Bags × 20 Tubes)	E1450-5050

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## Accessories

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	<a href="#">Rack 5.0ml tubes, blue</a>	5 Pcs. (1 Pack × 5 Pcs.)	E2325-5011
	<a href="#">Universal centrifuge adapter, for 5.0 ml tubes, for 15 ml rotor bores</a>	1 Pack × 8 Pcs.	E1450-9011
	<a href="#">Universal centrifuge adapter, for 5.0 ml tubes, for 50 ml rotor bores</a>	1 Pack × 8 Pcs.	E1450-9015
	<a href="#">Cap Lock for 5.0ml Centrifuge tubes</a>	100 Pcs. (10 Bags × 10 Pcs.)	E1450-9000
	<a href="#">Storage box 5.0ml tubes, blue</a>	3 Pcs. (1 Pack × 3 Pcs.)	E2325-5021

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)