



## 12-Place Rack for Upright Freezer, Pull Out Shelves for StarStore 100



### Product Highlights

- 12 Places (3 deep x 4 high)
- Racks with front access for use in upright freezers have drawers with handle
- For use with all standard storage boxes up to 140 x 155 x 55 mm
- Made from 85% stainless steel
- Made in a factory powered 100% by sustainable energy
- Article No. etched onto every rack for easy reordering

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	E6001-0134
Pack Size	1 Piece (1 Rack × 1 Piece)
Dimensions (W × D × H)	147 mm × 500 mm × 239 mm
Dimensions (W × D)	147 mm × 500 mm
External width	147 mm
Rack Type	Upright Freezer Rack
Material	Stainless Steel
Places	12 ( 3 deep x 4 high)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)









## More informations about 12-Place Rack for Upright Freezer, Pull Out Shelves for StarStore 100

Upright freezer racks made specifically for use with our popular StarStore 100 storage boxes (or any other storage box up to 140 x 155 x 55 mm). Keeping your freezer tidy and optimising storage space is an efficient use of energy. Racks with pull out shelves provide the convenience of quick access to your sample without the need to remove the rack from the freezer.



## All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	<a href="#">9-Place Rack for Upright Freezer, Pull Out Shelves for StarStore 100</a> <b>Places:</b> 9 ( 3 deep x 3 high)	1 Piece (1 Rack x 1 Piece)	E6001-0133
	 <a href="#">Create PDF Data Sheet</a>		
	<a href="#">12-Place Rack for Upright Freezer, Pull Out Shelves for StarStore 100</a> <b>Places:</b> 12 ( 3 deep x 4 high)	1 Piece (1 Rack x 1 Piece)	E6001-0134
	 <a href="#">Create PDF Data Sheet</a>		
	<a href="#">16-Place Rack for Upright Freezer, Pull Out Shelves for StarStore 100</a> <b>Places:</b> 16 ( 4 deep x 4 high)	1 Piece (1 Rack x 1 Piece)	E6001-0144
	 <a href="#">Create PDF Data Sheet</a>		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)