



Chest Freezer Racks for Standard Storage Boxes



Product Highlights

- Racks for use in chest freezers have a locking rod to keep boxes secure and a handle on top
- For use with all standard storage boxes up to 136 x 136 x 53 mm
- Laser cutting technology results in high quality units with no sharp edges
- Made from 85% stainless steel
- Made in a factory powered 100% by sustainable energy
- Article No. etched onto every rack for easy reordering

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	See variations
Rack Type	Chest Freezer Rack
Material	Stainless Steel

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Chest Freezer Racks for Standard Storage Boxes

Chest freezer racks for use with all standard storage boxes up to 136 x 136 x 53 mm. Keeping your freezer tidy and optimising storage space is an efficient use of energy. Racks for chest freezers feature a locking rod to keep your storage boxes secure and a handle on top for easy removal from the freezer

The following storage boxes are suitable for use with these racks:

- StarStore 81, StarStore 49, StarStore
- The 50 mm high cardboard storage boxes
- 81-place cryogenic storage boxes, 100-place cryogenic storage boxes
- StarRack 64



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 9-High Rack for Chest Freezer Places: 9 (9 high)  Create PDF Data Sheet	1 Piece (1 Rack × 1 Piece)	E6010-0009
 10-High Rack for Chest Freezer Places: 10 (10 high)  Create PDF Data Sheet	1 Piece (1 Rack × 1 Piece)	E6010-0010
 11-High Rack for Chest Freezer Places: 11 (11 high)  Create PDF Data Sheet	1 Piece (1 Rack × 1 Piece)	E6010-0011
 12-High Rack for Chest Freezer Places: 12 (12 high)  Create PDF Data Sheet	1 Piece (1 Rack × 1 Piece)	E6010-0012
 13-High Rack for Chest Freezer Places: 13 (13 high)  Create PDF Data Sheet	1 Piece (1 Rack × 1 Piece)	E6010-0013

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de