



## Screw Cap Microtubes with Natural Standard Caps (Sterile)



### Product Highlights

- Routinely tested to 20,000 x g
- Sterile and certified RNase, DNase, DNA, PCR Inhibitor and Endotoxin free.
- Tubes are supplied with pre-assembled natural standard caps
- Standard cap design allows simple and secure closure
- Tube made of polypropylene (PP)

## Product Applications

Suitable for many applications from sample storage in the gas phase of liquid nitrogen to freezer, through boiling to autoclave

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	See variations
Sterile	Yes
Autoclavable	No
Cap style	Natural standard cap
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
PCR inhibitor free	Yes
RNase free	Yes
Color	Natural
g-Force max.	20,000 x g

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)













## More informations about Screw Cap Microtubes with Natural Standard Caps (Sterile)

### Assembled sterile tubes & caps

A range of sterile and pre-assembled microtubes, available with a range of different cap styles. Supplied in convenient, resealable mini bags which stand up on your bench. Suitable for many applications from sample storage in the gas phase of liquid nitrogen through freezer to boiling.



## All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	<a href="#">0.5 ml Plain Skirted Tube, Natural Standard Screw Cap (Sterile)</a> <b>Volume:</b> 0.5 ml <b>Tube Base:</b> Skirted <b>Tube Style:</b> Plain   <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1405-2141
	<a href="#">1.5 ml Graduated Conical Tube, Natural Standard Screw Cap (Sterile)</a> <b>Volume:</b> 1.5 ml <b>Tube Base:</b> Conical <b>Tube Style:</b> Graduated   <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1415-2231
	<a href="#">1.5 ml Plain Skirted Tube, Natural Standard Screw Cap (Sterile)</a> <b>Volume:</b> 1.5 ml <b>Tube Base:</b> Skirted <b>Tube Style:</b> Plain   <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1415-2241
	<a href="#">2.0 ml Graduated Conical Tube, Natural Standard Screw Cap (Sterile)</a> <b>Volume:</b> 2.0 ml <b>Tube Base:</b> Conical <b>Tube Style:</b> Graduated   <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1420-2331
	<a href="#">2.0 ml Graduated Skirted Tube, Natural Standard Screw Cap (Sterile)</a> <b>Volume:</b> 2.0 ml <b>Tube Base:</b> Skirted <b>Tube Style:</b> Graduated   <a href="#">Create PDF Data Sheet</a>	250 Tubes (5 Bags × 50 Tubes)	E1420-2341

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)