



## 8-Strip PCR Tubes



### Product Highlights

- Manufactured from medical grade polypropylene
- Tubes can easily be split into smaller sections by snapping the strips apart
- Ultra thin walls ensure efficient thermal transfer and maximum yield
- PCR inhibitor free
- Certified RNase, DNase, DNA and Pyrogen-free
- Supplied in resealable bags

## Product Applications

PCR

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	See variations
Sterile	No
Autoclavable	Yes
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Human DNA free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Various
Suitable for PCR	yes
Suitable for Real-Time PCR (qPCR)	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)













## More informations about 8-Strip PCR Tubes

8-strip [PCR tubes](#). A range of different types of [PCR strip caps](#) are available separately, so you can mix and match tube and cap style to suit your needs.





## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">0.2 ml 8-Strip PCR Tubes, Natural</a> <b>Volume:</b> 0.2 ml <b>Color:</b> Natural <b>Pack Size:</b> 125 Tube Strips (1 Bag × 125 Tube Strips)  <a href="#">Create PDF Data Sheet</a>	125 Tube Strips (1 Bag × 125 Tube Strips)	I1402-3500
 <a href="#">0.2 ml 8-Strip PCR Tubes (Case Size), natural</a> <b>Volume:</b> 0.2 ml <b>Color:</b> Natural <b>Pack Size:</b> 1,250 Tube Strips (10 Bags × 125 Tube Strips)  <a href="#">Create PDF Data Sheet</a>	1,250 Tube Strips (10 Bags × 125 Tube Strips)	I1402-3500-C
 <a href="#">0.2 ml 8-Strip PCR Tubes, Mixed</a> <b>Volume:</b> 0.2 ml <b>Color:</b> Mixed <b>Pack Size:</b> 125 Tube Strips (5 Bags × 25 Tube Strips)  <a href="#">Create PDF Data Sheet</a>	125 Tube Strips (5 Bags × 25 Tube Strips)	I1402-3508
 <a href="#">0.2 ml 8-Strip PCR Tubes, White</a> <b>Volume:</b> 0.2 ml <b>Color:</b> White <b>Pack Size:</b> 125 Tube Strips (1 Bag × 125 Tube Strips)  <a href="#">Create PDF Data Sheet</a>	125 Tube Strips (1 Bag × 125 Tube Strips)	I1402-3509
 <a href="#">0.1 ml 8-Strip PCR Tubes, Low-Profile, Natural</a> <b>Volume:</b> 0.1 ml <b>Color:</b> Natural <b>Pack Size:</b> 125 Tube Strips (1 Bag × 125 Tube Strips)  <a href="#">Create PDF Data Sheet</a>	125 Tube Strips (1 Bag × 125 Tube Strips)	I1402-3600

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p><a href="#">0.1 ml 8-Strip PCR Tubes, Low-Profile, white</a> <b>Volume:</b> 0.1 ml <b>Color:</b> White <b>Pack Size:</b> 125 Tube Strips (1 Bag × 125 Tube Strips)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	125 Tube Strips (1 Bag × 125 Tube Strips)	I1402-3609

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)