



Cryogenic Lock Racks



Product Highlights

- 50-place polypropylene rack in a 5x10 array for skirted cryogenic vials
- Skirted vials lock into the rack for easy one-handed cap removal
- Anti-skid rubber feet
- Temperature use: -90 °C to 121 °C
- Autoclavable
- Dimensions (WxDxH): 200 x 100 x 25 mm

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	See variations
Pack Size	4 Pcs. (4 Pcs.)
Dimensions (W × D × H)	200 mm × 100 mm × 25 mm
Dimensions (W × D)	200 mm × 100 mm
External width	200 mm
Autoclavable	Yes
Material	Polypropylene
Color	Various
Places	50

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Cryogenic Lock Racks







50-place polypropylene rack in three different colors for skirted cryogenic vials. The vials lock into the rack for easy one-handed cap removal.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	50 Place Lock-Rack, Blue Color: Blue	4 Pcs. (1 Pack × 4 Pcs.)	E2350-4001
	 Create PDF Data Sheet		
	50 Place Lock-Rack, Red Color: Red	4 Pcs. (1 Pack × 4 Pcs.)	E2350-4004
	 Create PDF Data Sheet		
	50 Place Lock-Rack, Violet Color: Violet	4 Pcs. (1 Pack × 4 Pcs.)	E2350-4005
	 Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de