



12-Strip PCR Tubes



Product Highlights

- Manufactured from medical grade polypropylene
- Tubes can easily be split into smaller sections by snapping the strips apart
- Ultra thin walls ensure efficient thermal transfer and maximum yield
- PCR inhibitor free
- Certified RNase, DNase, DNA and Pyrogen-free
- Supplied in resealable bags

Product Applications

PCR

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	See variations
Sterile	No
Autoclavable	Yes
Volume	0.2 ml
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Human DNA free	Yes
PCR inhibitor free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Various
Suitable for PCR	yes
Suitable for Real-Time PCR (qPCR)	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 12-Strip PCR Tubes

8- and 12-strip PCR tubes. A range of different types of PCR strip caps are available separately, so you can mix and match tube and cap style to suit your needs.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 0.2 ml 12-Strip PCR Tubes, Natural Color: Natural Pack Size: 80 Tube Strips (1 Bag × 80 Tube Strips)  Create PDF Data Sheet	80 Tube Strips (1 Bag × 80 Tube Strips)	I1402-1200
 0.2 ml 12-Strip PCR Tubes, Natural (Case Size) Color: Natural Pack Size: 800 Tube Strips (10 Bags × 80 Tube Strips)  Create PDF Data Sheet	800 Tube Strips (10 Bags × 80 Tube Strips)	I1402-1200-C
 0.2 ml 12-Strip PCR Tubes, Mixed Color: Mixed Pack Size: 80 Tube Strips (5 Bags × 16 Tube Strips)  Create PDF Data Sheet	80 Tube Strips (5 Bags × 16 Tube Strips)	I1402-1208
 0.2 ml 12-Strip PCR Tubes Color: Blue Pack Size: 16 Tube Strips (1 Bag × 16 Tube Strips)  Create PDF Data Sheet	16 Tube Strips (1 Bag × 16 Tube Strips)	I1402-1201

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de