



## Semi-Skirted PCR Plates with Raised Rim



### Product Highlights

- Approximate maximum capacity for all 96-well plates is 350  $\mu\text{l}$  for standard height plates, or 200  $\mu\text{l}$  for low-profile plates
- Available as natural, white or black plates. White plates are ideal for qPCR, giving optimal signal-to-noise ratio for fluorescent based assays. Black plates help to minimize light diffusion.
- Specialist plates (for FAST<sup>®</sup> Systems, qPCR)

## Product Applications

PCR

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	See variations
Sterile	No
Autoclavable	Yes
Number of wells	96 wells
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
Free of heavy metals	Yes
PCR inhibitor free	Yes
RNase free	Yes
Frame Design	Semi-skirted
Suitable for PCR	yes
Suitable for Real-Time PCR (qPCR)	yes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about Semi-Skirted PCR Plates with Raised Rim









STARLAB Semi-Skirted PCR Plates feature a 7 mm high skirt around the complete plate edge.

The skirt improves the stability of the plate and makes them more suited e. g. for automated applications. Furthermore a labelling on the plate edge e. g. by barcode-labels is possible. Plates with elevated wells reduce the contamination risk by cross-contamination.

- > Autoclavable
- > Supplied in resealable bags
- > Convenient small pack sizes
- > Certified RNase, DNase, DNA and PCR Inhibitor free
- > Endotoxin free







## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p><a href="#">96-Well PCR Plate, Semi-Skirted, Raised Rim, Low Profile (for FAST® Systems), natural</a></p> <p><b>Profile:</b> Low (for FAST<sup>®</sup> Systems), raised rim  <b>Color:</b> Natural  <b>Pack Size:</b> 10 Plates (1 Box × 10 Plates)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	10 Plates (1 Box × 10 Plates)	E1403-7700
 <p><a href="#">96-Well PCR Plate, Semi-Skirted, Raised Rim, Low Profile (for FAST® Systems), natural (Case Size)</a></p> <p><b>Profile:</b> Low (for FAST<sup>®</sup> Systems), raised rim  <b>Color:</b> Natural  <b>Pack Size:</b> 100 Plates (10 Boxes × 10 Plates)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	100 Plates (10 Boxes × 10 Plates)	E1403-7700-C
 <p><a href="#">96-Well PCR Plate, Semi-Skirted, Raised Rim, Low Profile (for FAST® Systems), white</a></p> <p><b>Profile:</b> Low (for FAST<sup>®</sup> Systems), raised rim  <b>Color:</b> White  <b>Pack Size:</b> 10 Plates (1 Box × 10 Plates)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	10 Plates (1 Box × 10 Plates)	E1403-7709
 <p><a href="#">96-Well PCR Plate, Semi-Skirted, Raised Rim, Low Profile (for FAST® Systems), white (Case Size)</a></p> <p><b>Profile:</b> Low (for FAST<sup>®</sup> Systems), raised rim  <b>Color:</b> White  <b>Pack Size:</b> 100 Plates (10 Boxes × 10 Plates)</p> <p> <a href="#">Create PDF Data Sheet</a></p>	100 Plates (10 Boxes × 10 Plates)	E1403-7709-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
 Neuer Höltigbaum 38  
 22143 Hamburg  
 Email: [info@starlab.de](mailto:info@starlab.de)



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">96-Well PCR Plate, Semi-Skirted, Raised Rim, natural</a> <b>Profile:</b> Raised rim <b>Color:</b> Natural <b>Pack Size:</b> 10 Plates (1 Box × 10 Plates)  <a href="#">Create PDF Data Sheet</a>	10 Plates (1 Box × 10 Plates)	E1403-8200
 <a href="#">96-Well PCR Plate, Semi-Skirted, Raised Rim, natural (Case Size)</a> <b>Profile:</b> Raised rim <b>Color:</b> Natural <b>Pack Size:</b> 100 Plates (10 Boxes × 10 Plates)  <a href="#">Create PDF Data Sheet</a>	100 Plates (10 Boxes × 10 Plates)	E1403-8200-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)