



Semi-Skirted PCR Plates with Raised Rim



Product Highlights

- Approximate maximum capacity for all 96-well plates is 350 μ l for standard height plates, or 200 μ l for low-profile plates
- Available as natural, white or black plates. White plates are ideal for qPCR, giving optimal signal-to-noise ratio for fluorescent based assays. Black plates help to minimize light diffusion.
- Specialist plates (for FAST[®] Systems, qPCR)

Product Applications

PCR

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

| | |
|-----------------------------------|----------------|
| Art. No. | See variations |
| Sterile | No |
| Autoclavable | Yes |
| Number of wells | 96 wells |
| Material | Polypropylene |
| DNA free | Yes |
| DNase free | Yes |
| Free of endotoxins | Yes |
| Free of heavy metals | Yes |
| PCR inhibitor free | Yes |
| RNase free | Yes |
| Frame Design | Semi-skirted |
| Suitable for PCR | yes |
| Suitable for Real-Time PCR (qPCR) | yes |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Semi-Skirted PCR Plates with Raised Rim









STARLAB Semi-Skirted PCR Plates feature a 7 mm high skirt around the complete plate edge.

The skirt improves the stability of the plate and makes them more suited e. g. for automated applications. Furthermore a labelling on the plate edge e. g. by barcode-labels is possible. Plates with elevated wells reduce the contamination risk by cross-contamination.

- > Autoclavable
- > Supplied in resealable bags
- > Convenient small pack sizes
- > Certified RNase, DNase, DNA and PCR Inhibitor free
- > Endotoxin free







All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|-----------------------------------|--------------|
|  96-Well PCR Plate, Semi-Skirted, Raised Rim, Low Profile (for FAST® Systems), natural Profile: Low (for FAST® Systems), raised rim Color: Natural Pack Size: 10 Plates (1 Box × 10 Plates)  Create PDF Data Sheet | 10 Plates (1 Box × 10 Plates) | E1403-7700 |
|  96-Well PCR Plate, Semi-Skirted, Raised Rim, Low Profile (for FAST® Systems), natural (Case Size) Profile: Low (for FAST® Systems), raised rim Color: Natural Pack Size: 100 Plates (10 Boxes × 10 Plates)  Create PDF Data Sheet | 100 Plates (10 Boxes × 10 Plates) | E1403-7700-C |
|  96-Well PCR Plate, Semi-Skirted, Raised Rim, Low Profile (for FAST® Systems), white Profile: Low (for FAST® Systems), raised rim Color: White Pack Size: 10 Plates (1 Box × 10 Plates)  Create PDF Data Sheet | 10 Plates (1 Box × 10 Plates) | E1403-7709 |
|  96-Well PCR Plate, Semi-Skirted, Raised Rim, Low Profile (for FAST® Systems), white (Case Size) Profile: Low (for FAST® Systems), raised rim Color: White Pack Size: 100 Plates (10 Boxes × 10 Plates)  Create PDF Data Sheet | 100 Plates (10 Boxes × 10 Plates) | E1403-7709-C |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|-----------------------------------|--------------|
|  96-Well PCR Plate, Semi-Skirted, Raised Rim, natural Profile: Raised rim Color: Natural Pack Size: 10 Plates (1 Box × 10 Plates)  Create PDF Data Sheet | 10 Plates (1 Box × 10 Plates) | E1403-8200 |
|  96-Well PCR Plate, Semi-Skirted, Raised Rim, natural (Case Size) Profile: Raised rim Color: Natural Pack Size: 100 Plates (10 Boxes × 10 Plates)  Create PDF Data Sheet | 100 Plates (10 Boxes × 10 Plates) | E1403-8200-C |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de