



Piston grease for ErgoOne® E 8- and 12-channel pipettes 1250 µl



Product Highlights

- Reliable STARLAB quality.
- Supplied with battery, AC adapter, shelf holder (G9999-1009), reagent reservoir, 1 set of o-rings, 2 racks of tips and silicone grease

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

| | |
|--------------|---------------------------|
| Art. No. | G9999-1012 |
| Pack Size | 1 Piece (1 Box × 1 Piece) |
| For use with | ErgoOne® E Pipettes |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Piston grease for ErgoOne® E 8- and 12-channel pipettes 1250 µl

ErgoOne® E single and multi-channel pipettes provide reliable pipetting every day. ErgoOne® E are microprocessor controlled, battery-operated, piston-stroke pipettes, which have been optimised for ergonomic operation and ease of use.