



## Wall Mount for one ErgoOne® E pipette



### Product Highlights

- Reliable STARLAB quality.

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	G9999-1010
Pack Size	1 Piece (1 Box × 1 Piece)
For use with (SL)	ErgoOne® E Pipettes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about Wall Mount for one ErgoOne® E pipette

[ErgoOne](#)® E single and multi-channel [pipettes](#) provide reliable pipetting every day. [ErgoOne](#)® E are microprocessor controlled, battery-operated, piston-stroke [pipettes](#), which have been optimised for ergonomic operation and ease of use.