




Cardboard Storage Box for 50 ml Centrifuge Tubes, 16 places



Product Highlights

- Moisture-resistant coating
- 16-place box for 50 ml tubes
- 146 x 146 x 76 mm
- Removable dividers

Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	E2023-1650
Dimensions (W × D × H)	146 mm × 146 mm × 76 mm
Dimensions (W × D)	146 mm × 146 mm
Max. tube length	140 mm
Autoclavable	No
Pack Size	1 Piece (1 Box × 1 Piece)
Material	Coated Cardboard
Color	White
Places	16
Sustainability attributes	Recyclable
For use with (SL)	50 ml Centrifuge Tubes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Cardboard Storage Box for 50 ml Centrifuge Tubes, 16 places





White cardboard storage boxes with lift-off lid for the storage of 15 ml or 50 ml centrifuge tubes. Moisture-resistant coating. Removable dividers. Box dimensions (WxDxH): 146 x 146 x 76 mm.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 Cardboard Storage Box for 15 ml Centrifuge Tubes, 36 places Places: 36  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	E2023-1615
 Cardboard Storage Box for 50 ml Centrifuge Tubes, 16 places Places: 16  Create PDF Data Sheet	1 Piece (1 Box × 1 Piece)	E2023-1650

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de