




## 250 mm Cell Scraper with 10 mm Pivoting Blade



### Product Highlights

- Pivoting, 10 mm silicone blade
- Flexible ABS handle for easy handling in flasks
- Gamma sterilised (SAL 10<sup>-6</sup>)
- RNase, DNase and Pyrogen free
- Non cytotoxic and non-pyrogenic
- Individually wrapped and in convenient inner packs of 25

## Product Applications

Mechanical harvesting of cells in dishes, plates and flasks

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	CC7600-0250
Pack Size	100 Pcs. (1 Pack × 100 Pcs.)
Length	250 mm
Sterile	Yes
Blade Type	Pivoting
DNA free	Yes
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes
Color	Blue
Handle	Angled
Non-cytotoxic	Yes
For use with	small flasks

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about 250 mm Cell Scraper with 10 mm Pivoting Blade







### Rich harvest!

Cell scrapers comprise a gentle silicone blade and an ABS polymer handle. They are available in three sizes for cell harvesting in a variety of dishes, flasks or plates:

- › 11 mm fixed blade with 220 mm straight handle for plates and dishes
- › 10 mm pivoting blade with 250 mm single angle handle for small flasks
- › 12 mm pivoting blade with 320 mm dual angle handle for larger flasks



## All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	<a href="#">250 mm Cell Scraper with 10 mm Pivoting Blade</a> <b>Blade Type:</b> Pivoting <b>For use with:</b> small flasks   <a href="#">Create PDF Data Sheet</a>	100 Pcs. (1 Pack × 100 Pcs.)	CC7600-0250
	<a href="#">220 mm Cell Scraper with 11 mm Fixed Blade</a> <b>Blade Type:</b> Fixed <b>For use with:</b> plates, dishes   <a href="#">Create PDF Data Sheet</a>	100 Pcs. (1 Pack × 100 Pcs.)	CC7600-0220
	<a href="#">300 mm Cell Scraper with 20 mm Pivoting Blade (Sterile)</a> <b>Blade Type:</b> Pivoting <b>For use with:</b> large flasks   <a href="#">Create PDF Data Sheet</a>	150 Pcs. (1 Pack × 150 Pcs.)	CC7600-0300

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)