



## Cardboard Storage Boxes, 133 x 133 standard, 53 mm high



### Product Highlights

- Moisture resistant coating
- Standard 133x133 footprint, for use in most freezer racking units
- Removable dividers
- Made in Germany, in a factory powered by renewable energy

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	See variations
Dimensions (W × D × H)	133 mm × 133 mm × 53 mm
Dimensions (W × D)	133 mm × 133 mm
Autoclavable	No
Material	Coated Cardboard
Color	White
Sustainability attributes	Made in a factory powered by renewable energy. Made from 100% recycled cardboard. Recyclable after use.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about Cardboard Storage Boxes, 133 x 133 standard, 53 mm high

Starlab's Cardboard [Storage Boxes](#) are made from high-quality 100% recycled cardboard and coated with a water-repellent film. This combination offers safe and reliable storage at extremely low temperatures. The standard footprint of 133x133x53mm provides compatibility with most stainless steel [freezer racking](#). Ease your sample inventory with the writing areas on top of the lid, and on the side of the lid for visibility of [labelling](#) when stored in racking. Supplied pre-assembled, with removable grid insert. The boxes are fully recyclable after use.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)