




10 ml Finn/Socorex-Style Blue Tip (sterile)



Product Highlights

- From 1200 µl tips to 10 ml tips
- Made from polypropylene
- Compatible with a wide range of pipettes
- Different pack types available
- Certified
- DNase free, RNase free, DNA free, Pyrogen free

Product Applications

Pipetting of large volumes of liquids

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	I1052-0010
Sterile	Yes
Autoclavable	Yes
Filter tip	Non-filter
Volume	10 ml
Pack Size	250 Tips (10 Racks × 25 Tips)
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes
Graduations	No
Design Tips	Finn/Socorex-Style
For use with brand + product (SL)	Sartorius®: Proline Plus Socorex®: Acura 835F, Acura electro 935, Calibra digital Thermo Fisher Scientific Inc.: Finn Focus, Finnpipette Digital, Finnpipette F1/F2, Finnpipette Novus

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 10 ml Finn/Socorex-Style Blue Tip (sterile)











Large Volume Tips

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de





All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 10 ml Finn/Socorex-Style Blue Tip (non-sterile) Design Tips: Finn/Socorex-Style Sterile: No Pack Size: 100 Tips (1 Bag × 100 Tips)  Create PDF Data Sheet	100 Tips (1 Bag × 100 Tips)	I1052-0000
 10 ml Finn/Socorex-Style Blue Tip (non-sterile) Design Tips: Finn/Socorex-Style Sterile: No Pack Size: 250 Tips (10 Racks × 25 Tips)  Create PDF Data Sheet	250 Tips (10 Racks × 25 Tips)	I1052-0001
 10 ml Finn/Socorex-Style Blue Tip (sterile) Design Tips: Finn/Socorex-Style Sterile: Yes Pack Size: 250 Tips (10 Racks × 25 Tips)  Create PDF Data Sheet	250 Tips (10 Racks × 25 Tips)	I1052-0010
 10 ml Gilson-Style Graduated Macro Tip (non-sterile) Design Tips: Gilson-Style Sterile: No Pack Size: 100 Tips (1 Bag × 100 Tips)  Create PDF Data Sheet	100 Tips (1 Bag × 100 Tips)	I1054-0000
 10 ml Gilson-Style Graduated Macro Tip (non-sterile) Design Tips: Gilson-Style Sterile: No Pack Size: 250 Tips (10 Racks × 25 Tips)  Create PDF Data Sheet	250 Tips (10 Racks × 25 Tips)	I1054-0001

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p>10 ml Gilson-Style Graduated Macro Tip (sterile) Design Tips: Gilson-Style Sterile: Yes Pack Size: 250 Tips (10 Racks × 25 Tips)</p> <p> Create PDF Data Sheet</p>	250 Tips (10 Racks × 25 Tips)	I1054-0011

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de