




5 ml Eppendorf-Style Graduated Macro Tip (non-sterile)



Product Highlights


- From 1200 μ l tips to 10 ml tips
- Made from polypropylene
- Compatible with a wide range of pipettes
- Different pack types available
- Certified
- DNase free, RNase free, DNA free, Pyrogen free

Product Applications

Pipetting of large volumes of liquids

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

| | |
|-----------------------------------|--|
| Art. No. | I1053-0000 |
| Sterile | No |
| Autoclavable | Yes |
| Filter tip | Non-filter |
| Volume | 5 ml |
| Pack Size | 250 Tips (1 Bag × 250 Tips) |
| Material | Polypropylene |
| DNA free | Yes |
| DNase free | Yes |
| Pyrogen free | Yes |
| RNase free | Yes |
| Graduations | Yes |
| Design Tips | Eppendorf-Style |
| For use with brand + product (SL) | Abimed : Discovery Comfort, LABMATE Optima Biohit : eLine, ePet, mLine, Picus, Proline Eppendorf : Reference 2, Research plus, Research pro, Xplorer Gilson : Pipetman Concept, Pipetman P, Pipetman Ultra Hamilton : SofTouch HTL : Discovery Comfort, Labmate Pro, Optipette Nichiryo : 31N, NPLO, NPX Rainin : Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus Sartorius : eLine, ePet, mLine, Picus, Picus NxT, Proline, Proline Plus, Tacta Socorex : Acura 835F STARLAB : MicroOne, StarPet |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 5 ml Eppendorf-Style Graduated Macro Tip (non-sterile)













Large Volume Tips

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de













All Variations

| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|-------------------------------|------------|
|  5 ml Biohit-Style Tip (non-sterile) Design Tips: Biohit-Style Sterile: No Pack Size: 250 Tips (1 Bag × 250 Tips)  Create PDF Data Sheet | 250 Tips (1 Bag × 250 Tips) | E1009-7000 |
|  5 ml Biohit-Style-XL Tip (non-sterile) Design Tips: Biohit-Style-XL Sterile: No Pack Size: 250 Tips (1 Bag × 250 Tips)  Create PDF Data Sheet | 250 Tips (1 Bag × 250 Tips) | E1009-7001 |
|  5 ml Finn-Style Tip (non-sterile) Design Tips: Finn-Style Sterile: No Pack Size: 250 Tips (1 Bag × 250 Tips)  Create PDF Data Sheet | 250 Tips (1 Bag × 250 Tips) | I1009-5000 |
|  5 ml Gilson-Style Graduated Macro Tip (non-sterile) Design Tips: Gilson-Style Sterile: No Pack Size: 250 Tips (1 Bag × 250 Tips)  Create PDF Data Sheet | 250 Tips (1 Bag × 250 Tips) | I1050-0000 |
|  5 ml Gilson-Style Graduated Macro Tip (non-sterile) Design Tips: Gilson-Style Sterile: No Pack Size: 500 Tips (10 Racks × 50 Tips)  Create PDF Data Sheet | 500 Tips (10 Racks × 50 Tips) | I1050-0700 |
|  5 ml Gilson-Style Graduated Macro Tip (sterile) Design Tips: Gilson-Style Sterile: Yes Pack Size: 500 Tips (10 Racks × 50 Tips)  Create PDF Data Sheet | 500 Tips (10 Racks × 50 Tips) | I1050-0710 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



| PRODUCT NAME | PACKAGING SIZE | ART. NO. |
|--|--------------------------------|------------|
|  <p>5 ml Eppendorf-Style Graduated Macro Tip (non-sterile) Design Tips: Eppendorf-Style Sterile: No Pack Size: 250 Tips (1 Bag × 250 Tips)</p> <p> Create PDF Data Sheet</p> | 250 Tips (1 Bag × 250 Tips) | I1053-0000 |
|  <p>5 ml Eppendorf-Style Graduated Macro Tip (non-sterile) Design Tips: Eppendorf-Style Sterile: No Pack Size: 500 Tips (10 Racks × 50 Tips)</p> <p> Create PDF Data Sheet</p> | 500 Tips (10 Racks × 50 Tips) | I1053-0700 |
|  <p>5 ml Eppendorf-Style Graduated Macro Tip (sterile) Design Tips: Eppendorf-Style Sterile: Yes Pack Size: 500 Tips (10 Racks × 50 Tips)</p> <p> Create PDF Data Sheet</p> | 500 Tips (10 Racks × 50 Tips) | I1053-0710 |
|  <p>5.0 ml ErgoOne E/Finn-Style Pipette Tip (non-sterile), Loose Design Tips: Finn-Style Sterile: No Pack Size: 200 Tips (1 Bag × 200 Tips)</p> <p> Create PDF Data Sheet</p> | 200 Tips (1 Bag × 200 Tips) | E1070-2596 |
|  <p>5.0 ml ErgoOne E/Finn-Style Pipette Tip (non-sterile), Loose Design Tips: Finn-Style Sterile: No Pack Size: 1,000 Tips (5 Bags × 200 Tips)</p> <p> Create PDF Data Sheet</p> | 1,000 Tips (5 Bags × 200 Tips) | E1070-2601 |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
 Neuer Hölftigbaum 38
 22143 Hamburg
 Email: info@starlab.de