



Benchtop Support for ErgoOne® Pipettes



Product Highlights

- Easy storage of ErgoOne® pipettes on the lab bench
- Suitable for the use with ErgoOne® single-channel pipettes

Product Applications

Easy storage of ErgoOne® pipettes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	E3099-0002
Pack Size	5 Pcs. (1 Pack × 5 Pcs.)
For use with (SL)	ErgoOne® Pipettes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Benchtop Support for ErgoOne® Pipettes

The benchtop support allows single-channel pipettes to be temporarily stored on the bench between use, without risking contamination.