



## Benchtop Support for ErgoOne® Pipettes



### Product Highlights

- Easy storage of ErgoOne® pipettes on the lab bench
- Suitable for the use with ErgoOne® single-channel pipettes

## Product Applications

Easy storage of ErgoOne® pipettes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	E3099-0002
Pack Size	5 Pcs. (1 Pack × 5 Pcs.)
For use with (SL)	ErgoOne® Pipettes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about Benchtop Support for ErgoOne® Pipettes

The benchtop support allows single-channel [pipettes](#) to be temporarily stored on the bench between use, without risking contamination.