



ErgoOne® FAST Benchtop Support



Product Highlights

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

| | |
|-------------------|----------------------------------|
| Art. No. | P7166-6306 |
| Pack Size | 1 Piece (1 Bag x 1 Piece) |
| For use with (SL) | ErgoOne® FAST Pipette Controller |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about ErgoOne® FAST Benchtop Support

During development of the ErgoOne® FAST electronic pipette controller, a high pipetting speed was recognised as an important feature. The result is convincing! The ErgoOne® FAST can fill a 25 ml pipette in less than 4 seconds!

Every aspect of ErgoOne® FAST has been developed with ergonomics in mind. ErgoOne® FAST is one of the lightest pipette controllers on the market. It is designed for comfort and balances perfectly in the hand to provide fatigue-free pipetting, even during long periods.

This benchtop support fits perfectly with the ErgoOne® FAST.