



StarPet Filler (Manual)



Product Highlights

- Use with all plastic or glass pipettes, 1.0 ml to 25 ml
- Exceptionally smooth operation
- Autoclavable silicon pipette holder
- 1 year warranty

Product Applications

Transfer of liquid volumes from 1.0 to 25 ml

Serial dispensing of aliquots of different volumes

Resuspension of bacteria or cell pellets and nucleic acid precipitates

Aspiration of cell layer from, e.g., Ficoll® gradient

Phase extraction and removal of lower phases

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	E4866-2000
Volume range	25 - 25 ml
Volume	25 ml
Pack Size	1 Piece (1 Box × 1 Piece)
Operating mode	Manual

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about StarPet Filler (Manual)

Use with all plastic or glass pipettes, 1.0 ml to 25 ml. This manual pipette filler has an exceptionally smooth operation. Two dispense modes: one delivers fast filling with optimal accuracy for 10 ml to 25 ml pipettes; the other is for use with 1.0 ml to 10 ml pipettes where accuracy can sometimes be difficult to control. Compatible with Blow-Out, 'To Deliver' (TD) or gravity dispense pipettes.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de