



100 ml Serological Pipette (sterile), individually wrapped



Product Highlights

- Accuracy $\pm 2\%$
- Bidirectional graduations (excl. 2.0 ml) provide added applicability
- Non-cytotoxic and non-haemolytic USP Class VI polystyrene
- Gamma sterilized to SAL 10^{-6} (ISO11137)
- Certified RNase, DNase, latex and pyrogen free
- Colour-coded for easy identification
- Negative graduations (excl. 100 ml) allows additional working volume
- Supplied with a polyolefin filter plug
- Individually wrapped in traceable paper-plastic peel-back packs showing lot number and expiry date
- Specially designed carton with perforated tear-strip doubles as efficient, bench top dispenser
- ISO 9001 and ISO 13485 manufactured

Product Applications

Serial dispensing of aliquots of different volumes

Resuspension of bacteria or cell pellets and nucleic acid precipitates

Aspiration of cell layer from, e.g., Ficoll® gradient

Phase extraction and removal of lower phases

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	E4860-1100
Length	346.8 mm
Sterile	Yes
Volume range	1 - 100 ml
Volume	100 ml
Pack Size	50 Pcs. (1 Pack × 50 Pcs.)
Color code	Pink
Material	Polystyrene
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes
Graduations	Yes
Color	Clear
For use with (SL)	ErgoOne® FAST Pipette Controller

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 100 ml Serological Pipette (sterile), individually wrapped

Discover the possibilities!















STARLAB's new extended range of Serological [Pipettes](#) cover seven different volumes and, for the first time, includes a 100 ml pipette. The [pipettes](#) are constructed from polystyrene for sterile serological and [tissue culture](#) applications. Colour-coding on both the [pipettes](#) and the packaging provides easy identification. The black printed, easy to read ascending and descending graduations simplify fluid measuring and dispensing, and all [pipettes](#) (excl. 100 ml) feature negative graduations for extra working volume. Perfect for use with all popular pipette controllers, including [ErgoOne](#)[®] FAST.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 1 ml Serological Pipette (sterile), individually wrapped Volume: 1 ml  Create PDF Data Sheet	500 Pcs. (1 Pack × 500 Pcs.)	E4860-0001
 2 ml Serological Pipette (sterile), individually wrapped Volume: 2 ml  Create PDF Data Sheet	500 Pcs. (1 Pack × 500 Pcs.)	E4860-0002
 5 ml Serological Pipette (sterile), individually wrapped Volume: 5 ml  Create PDF Data Sheet	200 Pcs. (1 Pack × 200 Pcs.)	E4860-0005
 10 ml Serological Pipette (sterile), individually wrapped Volume: 10 ml  Create PDF Data Sheet	200 Pcs. (1 Pack × 200 Pcs.)	E4860-0010
 25 ml Serological Pipette (sterile), individually wrapped Volume: 25 ml  Create PDF Data Sheet	150 Pcs. (1 Pack × 150 Pcs.)	E4860-0025
 50 ml Serological Pipette (sterile), individually wrapped Volume: 50 ml  Create PDF Data Sheet	100 Pcs. (1 Pack × 100 Pcs.)	E4860-0050
 100 ml Serological Pipette (sterile), individually wrapped Volume: 100 ml  Create PDF Data Sheet	50 Pcs. (1 Pack × 50 Pcs.)	E4860-1100

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de