



Mini Dry Bath Metal Block for 12 x 1.5ml Tubes (Conical)



Product Highlights

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	N2400-4023
Pack Size	1 Piece (1 Box × 1 Piece)
Dimensions (W × D × H)	71 mm × 47 mm × 32 mm
Dimensions (W × D)	71 mm × 47 mm
Well diameter	10.8 mm
Weight	0.6 kg
Well depth	28.5 mm
System	Mini Dry Bath Incubator
For use with	12 x 1.5 ml microcentrifuge tubes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Mini Dry Bath Metal Block for 12 x 1.5ml Tubes (Conical)

Small, but powerful!

A small but powerful unit suitable for a range of applications. The digital microprocessor controls the preset temperature of the metal block and the fast ramping of temperature saves your valuable time. Supplied with a transparent lid for temperature uniformity, the unit can also be used as a water bath.

We offer a range of accessories below.