



## Mini Dry Bath Metal Block for 6 x 15ml Tubes



### Product Highlights

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	N2400-4025
Pack Size	1 Piece (1 Box × 1 Piece)
Dimensions (W × D × H)	71 mm × 47 mm × 75 mm
Dimensions (W × D)	71 mm × 47 mm
Well diameter	17.3 mm
Weight	0.6 kg
Well depth	70 mm
System	Mini Dry Bath Incubator
For use with (SL)	6 x 15 ml centrifuge tubes

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about Mini Dry Bath Metal Block for 6 x 15ml Tubes

### Small, but powerful!

A small but powerful unit suitable for a range of applications. The digital microprocessor controls the preset temperature of the metal block and the fast ramping of temperature saves your valuable time. Supplied with a transparent lid for temperature uniformity, the unit can also be used as a water bath.

We offer a range of accessories below.