



## Adapter Plate 96 x 0.2 ml (for use with Thermoblock S8012-0018)



### Product Highlights

- Adapter plate for use with Thermoblock S8012-0018
- Converts S8012-0018 for use with 0.2 ml Tubes or one PCR Plate

## Product Applications

Sample Preparation

Purification of plasmids, DNA, RNA

cDNA synthesis

Enzyme reactions (restriction digest, ligations, proteinase-K-digestion)

Denaturation of DNA, RNA, proteins

Labeling of DNA, RNA, proteins

Cultivation of bacteria in microcentrifuge tubes, centrifuge tubes and deepwell plates

Transformation of bacteria

Resuspension of pellets

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	S8012-0019
Pack Size	1 Piece (1 Box × 1 Piece)
Operating temperature	4 to 35 °C
Sample capacity	96 × 0.2 ml PCR tubes / 12 × 0.2 ml 8-Strip PCR tubes / 8 × 0.2 ml 12-Strip PCR tubes / 1 × 96 × 0.2 ml PCR plate
Speed Range	300 – 1,400 rpm; 800 rpm when loadings > 200 g

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about Adapter Plate 96 x 0.2 ml (for use with Thermoblock S8012-0018)

Adapter plate to convert Thermoblock S8012-0018 for use with 0.2 ml Tubes (single or strips) or one PCR Plate.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)