



## ISO-Rack for 24 x 1.5/2.0 ml Microcentrifuge Tubes



### Product Highlights

- ISO-Rack for 24 x 1.5/2.0 ml Tubes
- For use with Mixer HC Thermoblocks S8012-0011 and S8012-0012

## Product Applications

Sample Preparation

Purification of plasmids, DNA, RNA

cDNA synthesis

Enzyme reactions (restriction digest, ligations, proteinase-K-digestion)

Denaturation of DNA, RNA, proteins

Labeling of DNA, RNA, proteins

Cultivation of bacteria in microcentrifuge tubes, centrifuge tubes and deepwell plates

Transformation of bacteria

Resuspension of pellets

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

|                       |                                       |
|-----------------------|---------------------------------------|
| Art. No.              | S8012-0021                            |
| Pack Size             | 1 Piece (1 Box × 1 Piece)             |
| Operating temperature | 4 to 35 °C                            |
| Sample capacity       | 24 × 1.5/2.0 ml microcentrifuge tubes |
| Speed Range           | 300 – 1,400 rpm                       |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about ISO-Rack for 24 x 1.5/2.0 ml Microcentrifuge Tubes

ISO rack for use with 1.5/2.0ml tubes (Mixer HC Thermoblocks S8012-0011 and S8012-0012). The Iso-Rack is ideal for time-critical applications and high throughput sample preparation. Prepare your samples in your tubes, using the Iso-Rack as your bench rack. You can then transfer all 24 tubes to the thermoblock at the same time, in seconds!