



Protection filter with connective tubing



Product Highlights

Product Applications

Aspiration of Liquids
Cell-Culture
Sterile Filtration
De-gassing of buffers
Solidphase Extraction

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	N2400-9007
Pack Size	1 Piece (1 Box × 1 Piece)

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Protection filter with connective tubing

Protection Filter with Connective Tubing. Pore size of 0.2 μm provides protection from biological hazards. Can be autoclaved for up to 20 times.