



## StarStore 50 Jr., Natural



### Product Highlights

- 50-place polypropylene storage box with fixed 5x10 dividers for tubes up to 10mm
- Lift-off lid with moulded matrix
- Temperature use -90°C to 121°C
- Autoclavable
- Dimensions (WxDxH): 121 x 62 x 37 mm
- Max. tube height: 60 mm

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

|  |                          |
|--|--------------------------|
| Art. No.                                     | E2305-6009               |
| Pack Size                                    | 5 Pcs. (1 Pack × 5 Pcs.) |
| Dimensions (W × D × H)                       | 121 mm × 63 mm × 38 mm   |
| Dimensions (W × D)                           | 121 mm × 63 mm           |
| External width                               | 121 mm                   |
| Autoclavable                                 | Yes                      |
| Material                                     | Polypropylene            |
| Color  | Natural                  |
| Max. centrifuge tubes or screwcap microtubes | 50                       |
| Places                                       | 50                       |

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about StarStore 50 Jr., Natural

50-place polypropylene storage box with fixed, 5x10 dividers for tubes up to 10mm, eg. 0.5 ml microcentrifuge tubes or PCR tubes.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltingbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)