




## 100 µl Round Gel-Loading Filter Tip



### Product Highlights

- Flexible, long and narrow
- Loading samples onto gels (eg. polyacrylamide or agarose)
- Handling of small volumes
- Compatible with a wide range of pipettes
- Different pack types available
- Certified

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	E1022-0008
Sterile	No
Filter tip	Filter
Volume	100 µl
Pack Size	1,000 Tips (1 Bag × 1000 Tips)
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Pyrogen free	Yes
RNase free	Yes
Tip Type	round gel-loading

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about 100 µl Round Gel-Loading Filter Tip









All tips are certified RNase and DNase free. Packaging of these tips differs from the TipOne® range.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">10 µl Round Gel-Loading Filter Tip</a> <b>Volume:</b> 10 µl <b>Sterile:</b> Yes  <a href="#">Create PDF Data Sheet</a>	1,020 Tips (5 Racks × 204 Tips)	I1010-3800
 <a href="#">20 µl Round Gel-Loading Filter Tip (Sterile)</a> <b>Volume:</b> 20 µl <b>Sterile:</b> Yes  <a href="#">Create PDF Data Sheet</a>	1,020 Tips (5 Racks × 204 Tips)	E1021-0618
 <a href="#">100 µl Round Gel-Loading Filter Tip</a> <b>Volume:</b> 100 µl <b>Sterile:</b> No  <a href="#">Create PDF Data Sheet</a>	1,000 Tips (1 Bag × 1000 Tips)	E1022-0008
 <a href="#">100 µl Round Gel-Loading Filter Tip (Sterile)</a> <b>Volume:</b> 100 µl <b>Sterile:</b> Yes  <a href="#">Create PDF Data Sheet</a>	1,020 Tips (5 Racks × 204 Tips)	I1022-0810

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)